



ENDODONTICS



Endodontics

Endodontic Access	19
Rotary Instruments	24
Paper Points	26
Root Canal Sealer	27
Obturators	28
Warm vertical obturation solution	28
Gutta-Percha Points	28
Accessories	30

Access	Shaping and Finishing	Irrigation	Obturation

- More space for debridement & debris extraction
- Respect of the natural tooth anatomy
- Preservation of tooth integrity while allowing for appropriate irrigation, disinfection and obturation


For a truly smooth, controlled and efficient experience.

dentsplysirona.com/trunatomy

#trunatomy



TruNatomy™ Orifice Modifier



TruNatomy™ Glider



TruNatomy™ Prime Endodontic File



TruNatomy™ Prime Absorbent Points



TruNatomy™ Conform Fit Gutta-Percha Points



TruNatomy™ Irrigation Needle

**+
WE
KNOW
ENDO.**

© 2019 Dentsply Sirona Inc. Rx Only. BDR / B EN TNMY ADV 000 / 00 / 2019 – created 01/2019



Cavity X-Gates Drill

A unique instrument that combines Gates Glidden drills 1, 2, 3 and 4 in one single efficient instrument with maximum efficiency.



X-Gates RA 32mm

A000824000700

Gates Glidden Drill

Gates Glidden drills are used to enlarge the upper portion of the root canal. They can be used with any rotating contra-angle handpiece.



Size (mm)	1 ø0.5	2 ø0.7	3 ø0.9	4 ø1.1	5 ø1.3	6 ø1.5
Assorted Pack	x1	x1	x1	x1	x1	x1
Individual Pack	x6	x6	x6	x6	x6	x6

Gates Glidden Drill RA 28mm 1	A000823000112
Gates Glidden Drill RA 28mm 2	A000823000212
Gates Glidden Drill RA 28mm 3	A000823000312
Gates Glidden Drill RA 28mm 4	A000823000412
Gates Glidden Drill RA 28mm 5	A000823000512
Gates Glidden Drill RA 28mm 6	A000823000612
Gates Glidden Drill RA 28mm ASS.1-6	A000823090112
Gates Glidden Drill RA 32mm 1	A000824000112
Gates Glidden Drill RA 32mm 2	A000824000212
Gates Glidden Drill RA 32mm 3	A000824000312
Gates Glidden Drill RA 32mm 4	A000824000412
Gates Glidden Drill RA 32mm 5	A000824000512
Gates Glidden Drill RA 32mm 6	A000824000612
Gates Glidden Drill RA 32mm ASS.1-6	A000824090112

C+ Files

The ideal hand files for the initial instrumentation of root canals. C+ Files have stronger buckling resistance compared to K-Files, that provides easier location of the canal orifices and easier access to the apical third. There is also better negotiation of calcifications and a better tactile feel. C+ Files have a depth gauge making it easier to determine working length when working with X-rays. The pyramid shaped tip makes insertion easier during instrumentation and the square sectioned file allows for better resistance to distortion. In addition, the polished surface of the C+ File allows for better removal of debris as well as a smoother insertion of the instrument in the canal.



Tip ISO Size	006	008	010	015
Assorted Tips 008-015*	-	x2	x2	x2
Individual Tips	x6	x6	x6	x6

Length: 18 mm, 21 mm, 25 mm and 31 mm
*Not available in 18mm length

C+File ReadySteel 18mm 006	A 012X01800604
C+File ReadySteel 18mm 008	A 012X01800804
C+File ReadySteel 18mm 010	A 012X01801004
C+File ReadySteel 18mm 015	A 012X01801504
C+File ReadySteel 18mm 020	A 012X01802004
C+File ReadySteel 18mm Asstd	A 012X01890004
C+File ReadySteel 21mm 006	A 012X02100604
C+File ReadySteel 21mm 008	A 012X02100804
C+File ReadySteel 21mm 010	A 012X02101004
C+File ReadySteel 21mm 015	A 012X02101504
C+File ReadySteel 21mm 020	A 012X02102004
C+File ReadySteel 21mm Asstd	A 012X02190004
C+File ReadySteel 25mm 006	A 012X02500604
C+File ReadySteel 25mm 008	A 012X02500804
C+File ReadySteel 25mm 010	A 012X02501004
C+File ReadySteel 25mm 015	A 012X02501504
C+File ReadySteel 25mm 020	A 012X02502004
C+File ReadySteel 25mm Asstd	A 012X02590004



K-File

The Maillefer K-File has a good cutting efficiency for straight root canals. Non-slip handles provide excellent grip and tactile responsiveness from the canal.



Tip ISO Size	006	008	010	015	020
Assorted Tips 015-040	-	-	-	x1	x1
Individual Tips	x6	x6	x6	x6	x6

Tip ISO Size	025	030	035	040	045
Assorted Tips 015-040	x1	x1	x1	x1	-
Assorted Tips 045-080	-	-	-	-	x1
Individual Tips	x6	x6	x6	x6	x6

Tip ISO Size	050	055	060	070	080
Assorted Tips 045-080	x1	x1	x1	x1	x1
Individual Tips	x6	x6	x6	x6	x6

Tip ISO Size	090	100	110	120	130
Individual Tips	x6	x6	x6	x6	x6*

Tip ISO Size	140	Length:
Individual Tips	x6#	21 mm, 25 mm and 31 mm

K-File M-Access 21MM 008	A12MA02100812
K-File M-Access 21MM 010	A12MA02101012
K-File M-Access 21MM 015	A12MA02101512
K-File M-Access 21MM 020	A12MA02102012
K-File M-Access 21MM 025	A12MA02102512
K-File M-Access 21MM 030	A12MA02103012
K-File M-Access 21MM 035	A12MA02103512
K-File M-Access 21MM 040	A12MA02104012
K-File M-Access 21MM 045	A12MA02104512
K-File M-Access 21MM 050	A12MA02105012
K-File M-Access 21MM 055	A12MA02105512
K-File M-Access 21MM 060	A12MA02106012
K-File M-Access 21MM 070	A12MA02107012
K-File M-Access 21MM 080	A12MA02108012
K-File M-Access 21MM 090	A12MA02109012
K-File M-Access 21MM 100	A12MA02110012
K-File M-Access 21MM 110	A12MA02111012
K-File M-Access 21MM 120	A12MA02112012
K-File M-Access 21MM 130	A12MA02113012
K-File M-Access 21MM 140	A12MA02114012
K-File M-Access 25MM 140	A12MA02514012
K-File M-Access 25MM 15-40	A12MA02590012
K-File M-Access 25MM 45-80	A12MA02590112
K-File M-Access 25MM 90-140	A12MA02590212

K-File M-Access 31MM 006	A12MA03100612
K-File M-Access 31MM 008	A12MA03100812
K-File M-Access 31MM 010	A12MA03101012
K-File M-Access 31MM 015	A12MA03101512
K-File M-Access 31MM 020	A12MA03102012
K-File M-Access 31MM 025	A12MA03102512
K-File M-Access 31MM 030	A12MA03103012
K-File M-Access 31MM 035	A12MA03103512
K-File M-Access 31MM 040	A12MA03104012
K-File M-Access 31MM 045	A12MA03104512
K-File M-Access 31MM 050	A12MA03105012
K-File M-Access 31MM 055	A12MA03105512
K-File M-Access 31MM 060	A12MA03106012
K-File M-Access 31MM 070	A12MA03107012
K-File M-Access 31MM 080	A12MA03108012
K-File M-Access 31MM 090	A12MA03109012
K-File M-Access 31MM 100	A12MA03110012
K-File M-Access 31MM 110	A12MA03111012
K-File M-Access 31MM 120	A12MA03112012
K-File M-Access 31MM 130	A12MA03113012
K-File M-Access 31MM 140	A12MA03114012
K-File M-Access 31MM 15-40	A12MA03190012
K-File M-Access 31MM 45-80	A12MA03190112
K-File M-Access 31MM 90-140	A12MA03190212

Hedstroem File

The Hedstroem File offers aggressive cutting efficiency, but with less chance of file breakage because of the greater strength provided by its thicker core. Cuts only on the upward stroke, pulling debris away from the apex.



Tip ISO Size	008	010	015	020	025
Assorted Tips 015-040	-	-	x1	x1	x1
Individual Tips	x6	x6	x6	x6	x6

Tip ISO Size	030	035	040	045	050
Assorted Tips 015-040	x1	x1	x1	-	-
Assorted Tips 045-080	-	-	-	x1	x1
Individual Tips	x6	x6	x6	x6	x6

Hedstroem M-Access 21MM 010	A16MA02101012
Hedstroem M-Access 21MM 015	A16MA02101512
Hedstroem M-Access 21MM 020	A16MA02102012
Hedstroem M-Access 21MM 025	A16MA02102512
Hedstroem M-Access 21MM 030	A16MA02103012
Hedstroem M-Access 21MM 035	A16MA02103512
Hedstroem M-Access 21MM 040	A16MA02104012
Hedstroem M-Access 21MM 045	A16MA02104512
Hedstroem M-Access 21MM 050	A16MA02105012
Hedstroem M-Access 21MM 055	A16MA02105512
Hedstroem M-Access 21MM 060	A16MA02106012
Hedstroem M-Access 21MM 070	A16MA02107012
Hedstroem M-Access 21MM 080	A16MA02108012
Hedstroem M-Access 21MM 090	A16MA02109012
Hedstroem M-Access 21MM 100	A16MA02110012
Hedstroem M-Access 21MM 120	A16MA02112012
Hedstroem M-Access 21MM 130	A16MA02113012
Hedstroem M-Access 21MM 140	A16MA02114012
Hedstroem M-Access 21MM 15-40	A16MA02190012
Hedstroem M-Access 21MM 45-80	A16MA02190112
Hedstroem M-Access 21MM 90-140	A16MA02190212
Hedstroem M-Access 25MM 008	A16MA02500812
Hedstroem M-Access 25MM 010	A16MA02501012
Hedstroem M-Access 25MM 015	A16MA02501512
Hedstroem M-Access 25MM 020	A16MA02502012

Hedstroem M-Access 25MM 025	A16MA02502512
Hedstroem M-Access 25MM 030	A16MA02503012
Hedstroem M-Access 25MM 035	A16MA02503512
Hedstroem M-Access 25MM 040	A16MA02504012
Hedstroem M-Access 25MM 045	A16MA02504512
Hedstroem M-Access 25MM 050	A16MA02505012
Hedstroem M-Access 25MM 055	A16MA02505512
Hedstroem M-Access 25MM 060	A16MA02506012
Hedstroem M-Access 25MM 070	A16MA02507012
Hedstroem M-Access 25MM 080	A16MA02508012
Hedstroem M-Access 25MM 090	A16MA02509012
Hedstroem M-Access 25MM 100	A16MA02510012
Hedstroem M-Access 25MM 110	A16MA02511012
Hedstroem M-Access 25MM 120	A16MA02512012
Hedstroem M-Access 25MM 130	A16MA02513012
Hedstroem M-Access 25MM 140	A16MA02514012
Hedstroem M-Access 25MM 15-40	A16MA02590012
Hedstroem M-Access 25MM 45-80	A16MA02590112
Hedstroem M-Access 25MM 90-140	A16MA02590212
Hedstroem M-Access 31MM 008	A16MA03100812
Hedstroem M-Access 31MM 010	A16MA03101012
Hedstroem M-Access 31MM 015	A16MA03101512
Hedstroem M-Access 31MM 020	A16MA03102012
Hedstroem M-Access 31MM 025	A16MA03102512
Hedstroem M-Access 31MM 030	A16MA03103012
Hedstroem M-Access 31MM 035	A16MA03103512
Hedstroem M-Access 31MM 040	A16MA03104012
Hedstroem M-Access 31MM 045	A16MA03104512
Hedstroem M-Access 31MM 050	A16MA03105012
Hedstroem M-Access 31MM 055	A16MA03105512
Hedstroem M-Access 31MM 060	A16MA03106012
Hedstroem M-Access 31MM 070	A16MA03107012
Hedstroem M-Access 31MM 080	A16MA03108012
Hedstroem M-Access 31MM 090	A16MA03109012
Hedstroem M-Access 31MM 100	A16MA03110012
Hedstroem M-Access 31MM 110	A16MA03111012
Hedstroem M-Access 31MM 120	A16MA03112012
Hedstroem M-Access 31MM 130	A16MA03113012
Hedstroem M-Access 31MM 140	A16MA03114012
Hedstroem M-Access 31MM 15-40	A16MA03190012
Hedstroem M-Access 31MM 45-80	A16MA03190112
Hedstroem M-Access 31MM 90-140	A16MA03190212





Nitiflex® File

					
Tip ISO Size					
Assorted Tips	015	020	025	030	035
Individual Tips	x1	x1	x1	x1	x1
Individual Tips	x6	x6	x6	x6	x6

Tip ISO Size					
Assorted Tips	040	045	050	055	060
Individual Tips	x1	-	-	-	-
Individual Tips	x6	x6	x6	x6	x6

Length: 21 mm and 25 mm

Nitiflex File*Not CE*21mm 015	A012N02101512
Nitiflex File*Not CE*21mm 020	A012N02102012
Nitiflex File*Not CE*21mm 025	A012N02102512
Nitiflex File*Not CE*21mm 030	A012N02103012
Nitiflex File*Not CE*21mm 035	A012N02103512
Nitiflex File*Not CE*21mm 040	A012N02104012
Nitiflex File*Not CE*21mm 045	A012N02104512
Nitiflex File*Not CE*21mm 050	A012N02105012
Nitiflex File*Not CE*21mm 055	A012N02105512
Nitiflex File*Not CE*21mm 060	A012N02106012
Nitiflex File 21mm AS 015-040	A012N02190000
Nitiflex File*Not CE*25mm 015	A012N02501512
Nitiflex File*Not CE*25mm 020	A012N02502012
Nitiflex File*Not CE*25mm 025	A012N02502512
Nitiflex File*Not CE*25mm 030	A012N02503012
Nitiflex File*Not CE*25mm 035	A012N02503512
Nitiflex File*Not CE*25mm 040	A012N02504012
Nitiflex File*Not CE*25mm 045	A012N02504512
Nitiflex File*Not CE*25mm 050	A012N02505012
Nitiflex File*Not CE*25mm 055	A012N02505512
Nitiflex File*Not CE*25mm 060	A012N02506012
Nitiflex File 25mm AS 015-040	A012N02590000

FlexoFile®

These Maillefer stainless steel files have outstanding flexibility, cutting efficiency and a non-slip handle. Indicated for curved and narrow root canals, FlexoFiles have a unique non-cutting Batt™ Tip™ that avoids ledging and perforation of curved canals.



Tip ISO Size					
Assorted Tips 21, 25 & 31 mm 015-040	006	008	010	015	020
Individual Tips	-	-	-	x1	x1
Individual Tips	x6	x6	x6	x6	x6

Tip ISO Size					
Assorted Tips 21, 25 & 31 mm 015-040	025	030	035	040	045
Individual Tips	x1	x1	x1	x1	-
Assorted Tips 21, 25 & 31 mm 045-080	-	-	-	-	x1
Individual Tips	x6	x6	x6	x6	x6

Tip ISO Size					
Assorted Tips 21, 25 & 31 mm 045-080	050	055	060	070	080
Individual Tips	x1	x1	x1	x1	x1
Individual Tips	x6	x6	x6	x6	x6

Tip ISO Size						
Assorted Tips 31 mm 090-140	090	100	110	120	130	140
Individual Tips	x1	x1	x1	x1	x1	x1
Individual Tips	x6	x6	x6	x6	x6	x6

Length: 21 mm, 25 mm and 31 mm

K-Flexofile Readysteel Asstd 21mm 15-40	A012C02190004
K-Flexofile Readysteel Asstd 25mm 15-40	A012C02590004
K-Flexofile Readysteel Asstd 31mm 15-40	A012C03190004
K-FlexoFile Readysteel 21mm 015	A012C02101504
K-FlexoFile Readysteel 21mm 020	A012C02102004
K-FlexoFile Readysteel 21mm 025	A012C02102504
K-FlexoFile Readysteel 21mm 030	A012C02103004
K-FlexoFile Readysteel 21mm 035	A012C02103504
K-FlexoFile Readysteel 21mm 040	A012C02104004
K-FlexoFile Readysteel 25mm 015	A012C02501504
K-FlexoFile Readysteel 25mm 020	A012C02502004
K-FlexoFile Readysteel 25mm 025	A012C02502504
K-FlexoFile Readysteel 25mm 030	A012C02503004
K-FlexoFile Readysteel 25mm 035	A012C02503504
K-FlexoFile Readysteel 25mm 040	A012C02504004



Pro Taper Ultimate™ Hand use file

Specifically designed and heat-treated to optimize the role in the shaping procedure

For complicated cases or dentists who prefer manual treatment. Same properties as the rotaly ones sterile



B ST PULH 619 OSX	Pro Taper Ultimate™ Hand use SX 19 mm x6
B ST PULH 621 OSL	Pro Taper Ultimate™ Hand use Slider 21 mm x6
B ST PULH 625 OSL	Pro Taper Ultimate™ Hand use Slider 25 mm x6
B ST PULH 631 OSL	Pro Taper Ultimate™ Hand use Slider 31 mm x6
B ST PULH 621 OSH	Pro Taper Ultimate™ Hand use Shaper 21 mm x6
B ST PULH 625 OSH	Pro Taper Ultimate™ Hand use Shaper 25 mm x6
B ST PULH 631 OSH	Pro Taper Ultimate™ Hand use Shaper 31 mm x6
B ST PULH 621 OF1	Pro Taper Ultimate™ Hand use F1 21 mm x6
B ST PULH 625 OF1	Pro Taper Ultimate™ Hand use F1 25 mm x6
B ST PULH 631 OF1	Pro Taper Ultimate™ Hand use F1 31 mm x6
B ST PULH 621 OF2	Pro Taper Ultimate™ Hand use F2 21 mm x6
B ST PULH 625 OF2	Pro Taper Ultimate™ Hand use F2 25 mm x6
B ST PULH 631 OF2	Pro Taper Ultimate™ Hand use F2 31 mm x6
B ST PULH 621 OF3	Pro Taper Ultimate™ Hand use F3 21 mm x6
B ST PULH 625 OF3	Pro Taper Ultimate™ Hand use F3 25 mm x6
B ST PULH 631 OF3	Pro Taper Ultimate™ Hand use F3 31 mm x6
B ST PULH 621 OFX	Pro Taper Ultimate™ Hand use FX 21 mm x6
B ST PULH 625 OFX	Pro Taper Ultimate™ Hand use FX 25 mm x6
B ST PULH 631 OFX	Pro Taper Ultimate™ Hand use FX 31 mm x6
B ST PULH 621 FXL	Pro Taper Ultimate™ Hand use FXL 21 mm x6
B ST PULH 625 FXL	Pro Taper Ultimate™ Hand use FXL 25 mm x6
B ST PULH 631 FXL	Pro Taper Ultimate™ Hand use FXL 31 mm x6
B ST PULH 521 SQE	Pro Taper Ultimate™ Hand use Sequence 21 mm x6
B ST PULH 525 SQE	Pro Taper Ultimate™ Hand use Sequence 25 mm x6
B ST PULH 531 SQE	Pro Taper Ultimate™ Hand use Sequence 31 mm x6



Rotary Instruments

ProTaper Next® Treatment Files



ProTaper Next XA File Ster. 19mm 6-File	A080021900A03
ProTaper Next X1 File Ster. 21mm 3-File	A080322100103
ProTaper Next X1 File Ster. 25mm 3-File	A080322500103
ProTaper Next X1 File Ster. 31mm 3-File	A080323100103
ProTaper Next X2 File Ster. 21mm 3-File	A080322100203
ProTaper Next X2 File Ster. 25mm 3-File	A080322500203
ProTaper Next X2 File Ster. 31mm 3-File	A080323100203
ProTaper Next X3 File Ster. 21mm 3-File	A080322100303
ProTaper Next X3 File Ster. 25mm 3-File	A080322500303
ProTaper Next X3 File Ster. 31mm 3-File	A080323100303
ProTaper Next X4 File Ster. 21mm 3-File	A080322100403
ProTaper Next X4 File Ster. 25mm 3-File	A080322500403
ProTaper Next X4 File Ster. 31mm 3-File	A080323100403
ProTaper Next X5 File Ster. 21mm 3-File	A080322100503
ProTaper Next X5 File Ster. 25mm 3-File	A080322500503
ProTaper Next X5 File Ster. 31mm 3-File	A080323100503
ProTaper Next Asstd X1-X3 Ster.21 3-File	A08032219A003
ProTaper Next Asstd X1-X3 Ster.25 3-File	A08032259A003

ProTaper Next files have variable tapers for an optimised crown-down technique. Measurements shown indicate file length, tip size and the taper at 1 mm from the tip of the file.

Handles for hand instrumentation

6x Handles Assorted	A140900090000
---------------------	---------------

ProTaper® Gold

ProTaper Gold technology includes a series of “shaping” and “finishing files” that create the predictable ProTaper shapes, and a system proven and trusted by clinicians through the world. ProTaper Gold rotary files feature the same exact geometries as ProTaper Universal, but offer an increase in flexibility. This is especially important in the finishing files, when navigating challenging curves in the apical region.

File Size	19 mm	21 mm	25 mm	31 mm	
Shaping File X (SX)	●				x6
Shaping File 1 (S1)		●	●	●	x6
Shaping File 2 (S2)		●	●	●	x6
Finishing File 1 (F1) Tip size 20/ variable taper		●	●	●	x6
Finishing File 2 (F2) Tip size 25/ variable taper		●	●	●	x6
Finishing File 3 (F3) Tip size 30/ variable taper		●	●	●	x6
Finishing File 4 (F4) Tip size 40/ variable taper		●	●	●	x6
Finishing File 5 (F5) Tip size 50/ variable taper		●	●	●	x6

Blisters of 6 Files

ProTaper Gold Sx 19mm	A0410219G0103
ProTaper Gold S1 21mm	A0410221G0103
ProTaper Gold S1 25mm	A0410225G0103
ProTaper Gold S1 31mm	A0410231G0103
ProTaper Gold S2 21mm	A0410221G0203
ProTaper Gold S2 25mm	A0410225G0203
ProTaper Gold S2 31mm	A0410231G0203
ProTaper Gold F1 21mm	A0411221G0103
ProTaper Gold F1 25mm	A0411225G0103
ProTaper Gold F1 31mm	A0411231G0103
ProTaper Gold F2 21mm	A0411221G0203
ProTaper Gold F2 25mm	A0411225G0203
ProTaper Gold F2 31mm	A0411231G0203
ProTaper Gold F3 21mm	A0411221G0303
ProTaper Gold F3 25mm	A0411225G0303
ProTaper Gold F3 31mm	A0411231G0303
ProTaper Gold F4 21mm	A0411221G0403
ProTaper Gold F4 25mm	A0411225G0403
ProTaper Gold F4 31mm	A0411231G0403
ProTaper Gold F5 21mm	A0411221G0503

ProTaper Gold F5 25mm	A0411225G0503
ProTaper Gold F5 31mm	A0411231G0503
ProTaper Gold Asstd 21mm	A0409221G0103
ProTaper Gold Asstd 25mm	A0409225G0103
ProTaper Gold Asstd 31mm	A0409231G0103

TruNatomy®

An innovative, efficient and complete root canal treatment solution, that offers a superior performance, debridement and debris removal, with an increased respects the natural anatomy of the canal while preserving structural integrity.

TruNatomy Orifice Modifier 4-Files 16mm	BSTTNMM416OMD
TruNatomy Glider 4-Files 21mm	BSTTNMM421GLR
TruNatomy Glider 4-Files 25mm	BSTTNMM425GLR
TruNatomy Glider 4-Files 31mm	BSTTNMM431GLR
TruNatomy Small 4-Files 21mm	BSTTNMM421SML
TruNatomy Small 4-Files 25mm	BSTTNMM425SML
TruNatomy Small 4-Files 31mm	BSTTNMM431SML
TruNatomy Prime 4-Files 21mm	BSTTNMM421PRM
TruNatomy Prime 4-Files 25mm	BSTTNMM425PRM
TruNatomy Prime 4-Files 31mm	BSTTNMM431PRM
TruNatomy Medium 4-Files 21mm	BSTTNMM421MED
TruNatomy Medium 4-Files 25mm	BSTTNMM425MED
TruNatomy Medium 4-Files 31mm	BSTTNMM431MED
TruNatomy Sequence 3-Files 21mm MU	BSTTNMM321SQC
TruNatomy Sequence 3-Files 25mm MU	BSTTNMM325SQC
TruNatomy Sequence 3-Files 31mm MU	BSTTNMM331SQC
TruNatomy Assort 5-Files 16-21mm MU	BSTTNMM521AST
TruNatomy Assort 5-Files 16-25mm MU	BSTTNMM525AST
TruNatomy Assort 5-Files 16-31mm MU	BSTTNMM531AST

Pro Taper Ultimate™

Specifically designed and heat-treated to optimize the role in the shaping procedure



B ST PULR 619 OSX	Pro Taper Ultimate™ SX 19 mm x6
B ST PULR 321 OSL	Pro Taper Ultimate™ Slider 21 mm x3
B ST PULR 325 OSL	Pro Taper Ultimate™ Slider 25 mm x3
B ST PULR 331 OSL	Pro Taper Ultimate™ Slider 31 mm x3
B ST PULR 621 OSL	Pro Taper Ultimate™ Slider 21 mm x6
B ST PULR 625 OSL	Pro Taper Ultimate™ Slider 25 mm x6
B ST PULR 631 OSL	Pro Taper Ultimate™ Slider 31 mm x6
B ST PULR 321 OSH	Pro Taper Ultimate™ Shaper 21 mm x3
B ST PULR 325 OSH	Pro Taper Ultimate™ Shaper 25 mm x3
B ST PULR 331 OSH	Pro Taper Ultimate™ Shaper 31 mm x3
B ST PULR 621 OSH	Pro Taper Ultimate™ Shaper 21 mm x6
B ST PULR 625 OSH	Pro Taper Ultimate™ Shaper 25 mm x6
B ST PULR 631 OSH	Pro Taper Ultimate™ Shaper 31 mm x6
B ST PULR 321 OF1	Pro Taper Ultimate™ F1 21 mm x3
B ST PULR 325 OF1	Pro Taper Ultimate™ F1 25 mm x3
B ST PULR 331 OF1	Pro Taper Ultimate™ F1 31 mm x3
B ST PULR 621 OF1	Pro Taper Ultimate™ F1 21 mm x6
B ST PULR 625 OF1	Pro Taper Ultimate™ F1 25 mm x6
B ST PULR 631 OF1	Pro Taper Ultimate™ F1 31 mm x6
B ST PULR 321 OF2	Pro Taper Ultimate™ F2 21 mm x3
B ST PULR 325 OF2	Pro Taper Ultimate™ F2 25 mm x3
B ST PULR 331 OF2	Pro Taper Ultimate™ F2 31 mm x3
B ST PULR 621 OF2	Pro Taper Ultimate™ F2 21 mm x6
B ST PULR 625 OF2	Pro Taper Ultimate™ F2 25 mm x6
B ST PULR 631 OF2	Pro Taper Ultimate™ F2 31 mm x6
B ST PULR 321 OF3	Pro Taper Ultimate™ F3 21 mm x3
B ST PULR 325 OF3	Pro Taper Ultimate™ F3 25 mm x3
B ST PULR 331 OF3	Pro Taper Ultimate™ F3 31 mm x3
B ST PULR 621 OF3	Pro Taper Ultimate™ F3 21 mm x6
B ST PULR 625 OF3	Pro Taper Ultimate™ F3 25 mm x6
B ST PULR 631 OF3	Pro Taper Ultimate™ F3 31 mm x6
B ST PULR 621 OFX	Pro Taper Ultimate™ FX 21 mm x6
B ST PULR 625 OFX	Pro Taper Ultimate™ FX 25 mm x6
B ST PULR 631 OFX	Pro Taper Ultimate™ FX 31 mm x6
B ST PULR 621 FXL	Pro Taper Ultimate™ FXL 21 mm x6
B ST PULR 625 FXL	Pro Taper Ultimate™ FXL 25 mm x6
B ST PULR 631 FXL	Pro Taper Ultimate™ FXL 31 mm x6
B ST PULR 521 SQE	Pro Taper Ultimate™ Sequence 21 mm x6
B ST PULR 525 SQE	Pro Taper Ultimate™ Sequence 25 mm x6
B ST PULR 531 SQE	Pro Taper Ultimate™ Sequence 31 mm x6



ProRinse®

Irrigation probe with a side port dispersal to improve debris removal



ProRinse®	A 0821
PRORINSE 23GA 100PK	A082100023100
PRORINSE 24GA 100PK	A082100024100
PRORINSE 25GA 100PK	A082100025100
PRORINSE 28GA 100PK	A082100028100
PRORINSE 30GA 100PK	A082100030100



Trunatomy® Irrigation Needle	BNSTNIRD2230G
• TruNatomy® Irrigation Needle 30 Gauge x40	

Glyde™

Glyde File Prep Gel is used to facilitate the cleaning and shaping of the root canal by lifting debris coronally, while it simultaneously lubricates the file. Glyde File Prep Gel effervesces when it reacts with sodium hypochlorite. This 'bubbling' action lifts dentinal mud and necrotic tissue for easy removal. In addition, the EDTA in Glyde assists to remove bacteria and helps to leave the dentinal tubes and lateral canals ready for obturation.



Glyde Syringe Kit (3 x 3 mL Syringes)	A090100000000
---------------------------------------	---------------

Paper Points

Branded Paper Points, specifically match the shape of matching branded files. They are colour coded for easy identification and come sterile.

ProTaper Next® Paper Points

ProTaper Next Paper Points (180 pack)

X2		A022Q00000200
X3		A022Q00000300
X2 & X3 Asstd (90 of each)		A022Q00090000
X4 & X5 Asstd (90 of each)		A022Q00090100

ProTaper Gold® Paper Points

ProTaper Gold Paper Points (180 pack)

F1		A241W00000103
F2		A241W00000203
F3		A241W00000303
F4		A241W00000403
F5		A241W00000503
F1-5		A241W00090103

WaveOne® Gold Absorbent Points

Sterile Paper Points (Box 180)

Small		A175W00000S03
Primary		A175W00000P03
Medium		A175W00000M03
Large		A175W00000L03

Trunatomy® Absorbent Points

Trunatomy® Absorbent Points (180 pack)

Small		BSTTNPP000SML
Prime		BSTTNPP000PRM
Medium		BSTTNPP000MED
Assortment		BSTTNPP000AST

ProTaper Ultimate™ Absorbent Points

F1		B ST PULP MPP OF1
F2		B ST PULP MPP OF2
F3		B ST PULP MPP OF3
FX		B ST PULP MPP OFX
FXL		B ST PULP MPP FXL
F1-F3		B ST PULP MPP AST

Root Canal Sealer

AH Plus®

AH Plus is a two component paste root canal sealer, based on epoxide-amine resins. This easy to mix sealer conforms closely to the walls of the prepared root canal and provides long-term dimensional stability with minimal shrinkage upon setting. AH Plus provides good biocompatibility, radiopacity, colour stability and ease of removal. Handling is quick and easy. Available in paste to paste mix or Automix Syringe (AH Plus Jet).



AH Plus Intro Package	60620112
-----------------------	----------

- 1 Tube of Paste A, (4 mL)
- 1 Tube of Paste B, (4 mL)



AH Plus Bioceramic Sealer Starter Kit	BKTAHBS 000 000
Syringe Refill	B00AHBS 100 03G
Cannula Refill	B00AHBS 020 24G

Proroot MTA



Proroot MTA 4X0.5G	A040500000300
--------------------	---------------



Warm vertical obturation solution

Gutta-Smart™ Cordless Obturation Device

Fill canals conveniently and reliably with Gutta-Smart cordless obturation device and Conform Fit gutta-percha master cones. Designed to work together, they feature the same low-temperature flow characteristics and radiopacity for a total warm vertical obturation solution.



Gutta-Smart™ Kit Contents

B00GSIQ000KIT

- 1 Flow Handpiece
- 1 Pack Handpiece
- 1 Base Station
- 2 Gutta-Percha Cartridges
- 3 Electric Heat Pluggers
- 1 Heating Element
- 1 Bending Tool
- 1 Power Adapter + 4 Universal Plugs

Gutta-Smart™ Conform Fit™ Gutta-Percha Cartridges



20 gauge 10-pack	BNSGSIQA0020G
23 gauge 10-pack	BNSGSIQA0023G
25 gauge 10-pack	BNSGSIQA0025G

Gutta-Smart™ Thermal Response Tip



1 TIP

B00GSIQ100000

Gutta-Smart™ Electric Heat Plugger



Small (40/.025), 1 each	B00GSIQ100SML
Medium (50/.05), 1 each	B00GSIQ100MED
Large (60/.06), 1 each	B00GSIQ100LRG

Gutta-Percha Points

Conform Fit™ Gutta-Percha

ProTaper Next™, ProTaper® Gold and TruNatomy™ Conform Fit™ Gutta-Percha WaveOne™ match the shape of their systems' rotary files and are colour coded for easy identification.

ProTaper Next® Conform Fit™ Gutta-Percha

ProTaper Next® Conform Fit™ Gutta-Percha (60 pack)



Red	X2	60-pack	B00PNGPF000X2
Blue	X3	60-pack	B00PNGPF000X3
Black	X4	60-pack	B00PNGPF00X23
Yellow	X5	60-pack	B00PNGPF00X45

ProTaper Gold® Conform Fit™ Gutta-Percha


ProTaper Gold® Conform Fit™ Gutta-Percha (60 pack)



Yellow	F1	60-pack	B00PGGPF000F1
Red	F2	60-pack	B00PGGPF000F2
Blue	F3	60-pack	B00PGGPF000F3
Assorted	F2/F3	60-pack	B00PGGPF00F23
Assorted	F4/F5	60-pack	B00PGGPF00F45

WaveOne® Gold Conform Fit™ Gutta-Percha

WaveOne® Gold Conform Fit™ Gutta-Percha (60 pack)



Yellow	Small	60-pack	B00WGGPF00SML
Red	Primary	60-pack	B00WGGPF00PRI
Green	Medium	60-pack	B00WGGPF00MED
White	Large	60-pack	B00WGGPF00LRG
Assorted	S/P/M/L	60-pack	B00WGGPF00AST

TruNatomy® Conform Fit™ Gutta-Percha



Yellow	Small	60-pack	B00TNGPF00SML
Red	Prime	60-pack	B00TNGPF00PRM
Green	Medium	60-pack	B00TNGPF00MED
Assorted	S/P/M	60-pack	B00TNGPF00AST

Pro Taper Ultimate™ Confit™ Gutta-Purcha (60 pack)

Yellow	F1	60-pack	B NS PULG FGP OF1
Red	F2	60-pack	B NS PULG FGP OF2
Blue	F3	60-pack	B NS PULG FGP OF3
Green	FX	60-pack	B NS PULG FGP OFX
Yellow	FXL	60-pack	B NS PULG FGP FXL
Assorted	F1-F3	60-pack	B NS PULG FGP AST

Gutta-Percha Points

Dentsply Sirona Gutta-Percha Points are dense, highly radio-opaque and flexible. They can be used with any condensation technique to fill the root canal.

Gutta Percha Auxiliary XF	A022L00000200
Gutta Percha Auxiliary F	A022L00000300
Gutta Percha Auxiliary M	A022L00000400


Gutta-Percha Points – ISO Sizes

GP-Points Pink 015	A022C00001500
GP-Points Pink 020	A022C00002000
GP-Points Pink 025	A022C00002500
GP-Points Pink 030	A022C00003000
GP-Points Pink 035	A022C00003500
GP-Points Pink 040	A022C00004000
GP-Points Pink 045	A022C00004500
GP-Points Pink 050	A022C00005000
GP-Points Pink 055	A022C00005500
GP-Points Pink 060	A022C00006000
GP-Points Pink 070	A022C00007000
GP-Points Pink 080	A022C00008000






Accessories



Endo-M-Bloc
Dr. P. Guldener

A018400000000



Mini-Endo-Bloc

A032700000000



Clean-Stand

A268P00000000

Clean-Stand Sponges

A268S00000000

MAILLEFER

ProTaper Ultimate™


Performance unlimited



SHAPING



CLEANING



OBTURATION

Dentsply
Sirona

It is all about passion.

Endodontics expertise requires performance.

Introducing a new solution combining the latest generation of ProTaper files, enhanced disinfection and dedicated obturation that works seamlessly together.

Shaping: 1 Slider-Shaper-Finishers sequence to cover a full range of anatomical situations.

Cleaning: More than 25% better cleaning efficacy vs. sodium hypochlorite soak without activation.*

Obturation: Obturation with ideal fit 100% of the time.

dentsplysirona.com/protaperultimate

#endopassion

*Only applicable for EndoActivator®

PART OF A

