



Dental Lab

Cercon
Lucitone

56
63





Cercon®
Zirconium Oxide

Makes you smile
Information for the dental laboratory



Cercon XT 98mm zirconia disk

Cercon extra translucent zirconia for enhanced aesthetics



Cercon xt A1 disk 98 12	5366111012
Cercon xt A1 disk 98 14	5366111014
Cercon xt A1 disk 98 18	5366111018
Cercon xt A1 disk 98 25	5366111025
Cercon xt A2 disk 98 12	5366111112
Cercon xt A2 disk 98 14	5366111114
Cercon xt A2 disk 98 18	5366111118
Cercon xt A2 disk 98 25	5366111125
Cercon xt A3 disk 98 12	5366111212
Cercon xt A3 disk 98 14	5366111214
Cercon xt A3 disk 98 18	5366111218
Cercon xt A3 disk 98 25	5366111225
Cercon xt A3.5 disk 98 12	5366111312
Cercon xt A3.5 disk 98 14	5366111314
Cercon xt A3.5 disk 98 18	5366111318
Cercon xt A3.5 disk 98 25	5366111325
Cercon xt A4 disk 98 12	5366111412
Cercon xt A4 disk 98 14	5366111414
Cercon xt A4 disk 98 18	5366111418
Cercon xt A4 disk 98 25	5366111425
Cercon xt B1 disk 98 12	5366111512
Cercon xt B1 disk 98 14	5366111514
Cercon xt B1 disk 98 18	5366111518
Cercon xt B1 disk 98 25	5366111525
Cercon xt B2 disk 98 12	5366111612
Cercon xt B2 disk 98 14	5366111614
Cercon xt B2 disk 98 18	5366111618
Cercon xt B2 disk 98 25	5366111625
Cercon xt B3 disk 98 12	5366111712
Cercon xt B3 disk 98 14	5366111714
Cercon xt B3 disk 98 18	5366111718
Cercon xt B3 disk 98 25	5366111725
Cercon xt B4 disk 98 12	5366111812
Cercon xt B4 disk 98 14	5366111814
Cercon xt B4 disk 98 18	5366111818
Cercon xt B4 disk 98 25	5366111825

Cercon xt C1 disk 98 12	5366111912
Cercon xt C1 disk 98 14	5366111914
Cercon xt C1 disk 98 18	5366111918
Cercon xt C1 disk 98 25	5366111925
Cercon xt C2 disk 98 12	5366112012
Cercon xt C2 disk 98 14	5366112014
Cercon xt C2 disk 98 18	5366112018
Cercon xt C2 disk 98 25	5366112025
Cercon xt C3 disk 98 12	5366112112
Cercon xt C3 disk 98 14	5366112114
Cercon xt C3 disk 98 18	5366112118
Cercon xt C3 disk 98 25	5366112125
Cercon xt C4 disk 98 12	5366112212
Cercon xt C4 disk 98 14	5366112214
Cercon xt C4 disk 98 18	5366112218
Cercon xt C4 disk 98 25	5366112225
Cercon xt D2 disk 98 12	5366112312
Cercon xt D2 disk 98 14	5366112314
Cercon xt D2 disk 98 18	5366112318
Cercon xt D2 disk 98 25	5366112325
Cercon xt D3 disk 98 12	5366112412
Cercon xt D3 disk 98 14	5366112414
Cercon xt D3 disk 98 18	5366112418
Cercon xt D3 disk 98 25	5366112425
Cercon xt D4 disk 98 12	5366112512
Cercon xt D4 disk 98 14	5366112514
Cercon xt D4 disk 98 18	5366112518
Cercon xt D4 disk 98 25	5366112525
Cercon xt White disk 98 12	5366112612
Cercon xt White disk 98 14	5366112614
Cercon xt White disk 98 18	5366112618
Cercon xt White disk 98 25	5366112625
Cercon xt BL2 disk 98 12	5366112712
Cercon xt BL2 disk 98 14	5366112714
Cercon xt BL2 disk 98 18	5366112718
Cercon xt BL2 disk 98 25	5366112725



Cercon® ht – high translucent zirconia

The new zirconia standard for color accuracy: True Color Technology

Cercon ht with True Color Technology sets new standards for zirconia when it comes to reproducing the classic 16 VITA® shades. Cercon ht provides the highest level of color accuracy based on our decades of know-how in ceramic shade-matching, our specially selected pigments and our own mixing and shade formulations – and that for almost all popular “open” milling units.

dentsplysirona.com



* VITA is a trademark of VITA Zahnfabrik

Cercon HT 98mm zirconia disk

Cercon extra translucent zirconia for enhanced aesthetics



Cercon ht A1 disk 98 12	5366091012
Cercon ht A1 disk 98 14	5366091014
Cercon ht A1 disk 98 18	5366091018
Cercon ht A1 disk 98 25	5366091025
Cercon ht A2 disk 98 12	5366091112
Cercon ht A2 disk 98 14	5366091114
Cercon ht A2 disk 98 18	5366091118
Cercon ht A2 disk 98 25	5366091125
Cercon ht A3 disk 98 12	5366091212
Cercon ht A3 disk 98 14	5366091214
Cercon ht A3 disk 98 18	5366091218
Cercon ht A3 disk 98 25	5366091225
Cercon ht A3.5 disk 98 12	5366091312
Cercon ht A3.5 disk 98 14	5366091314
Cercon ht A3.5 disk 98 18	5366091318
Cercon ht A3.5 disk 98 25	5366091325
Cercon ht A4 disk 98 12	5366091412
Cercon ht A4 disk 98 14	5366091414
Cercon ht A4 disk 98 18	5366091418
Cercon ht A4 disk 98 25	5366091425
Cercon ht B1 disk 98 12	5366091512
Cercon ht B1 disk 98 14	5366091514
Cercon ht B1 disk 98 18	5366091518
Cercon ht B1 disk 98 25	5366091525
Cercon ht B2 disk 98 12	5366091612
Cercon ht B2 disk 98 14	5366091614
Cercon ht B2 disk 98 18	5366091618
Cercon ht B2 disk 98 25	5366091625
Cercon ht B3 disk 98 12	5366091712
Cercon ht B3 disk 98 14	5366091714
Cercon ht B3 disk 98 18	5366091718
Cercon ht B3 disk 98 25	5366091725
Cercon ht B4 disk 98 12	5366091812
Cercon ht B4 disk 98 14	5366091814
Cercon ht B4 disk 98 18	5366091818
Cercon ht B4 disk 98 25	5366091825

Cercon ht C1 disk 98 12	5366091912
Cercon ht C1 disk 98 14	5366091914
Cercon ht C1 disk 98 18	5366091918
Cercon ht C1 disk 98 25	5366091925
Cercon ht C2 disk 98 12	5366092012
Cercon ht C2 disk 98 14	5366092014
Cercon ht C2 disk 98 18	5366092018
Cercon ht C2 disk 98 25	5366092025
Cercon ht C3 disk 98 12	5366092112
Cercon ht C3 disk 98 14	5366092114
Cercon ht C3 disk 98 18	5366092118
Cercon ht C3 disk 98 25	5366092125
Cercon ht C4 disk 98 12	5366092212
Cercon ht C4 disk 98 14	5366092214
Cercon ht C4 disk 98 18	5366092218
Cercon ht C4 disk 98 25	5366092225
Cercon ht D2 disk 98 12	5366092312
Cercon ht D2 disk 98 14	5366092314
Cercon ht D2 disk 98 18	5366092318
Cercon ht D2 disk 98 25	5366092325
Cercon ht D3 disk 98 12	5366092412
Cercon ht D3 disk 98 14	5366092414
Cercon ht D3 disk 98 18	5366092418
Cercon ht D3 disk 98 25	5366092425
Cercon ht D4 disk 98 12	5366092512
Cercon ht D4 disk 98 14	5366092514
Cercon ht D4 disk 98 18	5366092518
Cercon ht D4 disk 98 25	5366092525
Cercon ht BL disk 98 12	5366092612
Cercon ht BL disk 98 14	5366092614
Cercon ht BL disk 98 18	5366092618
Cercon ht BL disk 98 25	5366092625
Cercon ht BL2 disk 98 12	5366092712
Cercon ht BL2 disk 98 14	5366092714
Cercon ht BL2 disk 98 18	5366092718
Cercon ht BL2 disk 98 25	5366092725



Cercon® xt ML – Extra translucent multilayer zirconia

Layers ahead of the rest

Take monolithic restorations to the next level

- **Perfect Lifelike Gradient** – Matching natural tooth aesthetics without the need to veneer or stain, especially for the anterior region (flexural strength: 750 MPa)
- **Guaranteed Color Accuracy*** – Unparalleled color matching covering all 16 VITA** shades and BL2
- **Easy to Use** – Easy to nest and same milling and firing protocols for all Cercon zirconia discs

dentsplysirona.com



* Based on our trusted True Color Technology
** VITA is a trademark of VITA Zahnfabrik

Cercon® xt ML– Disc 98 mm

Cercon extra translucent zirconia for enhanced aesthetics



Cercon xt ML A1 disk 98 14	5366115014
Cercon xt ML A1 disk 98 18	5366115018
Cercon xt ML A1 disk 98 25	5366115025
Cercon xt ML A2 disk 98 14	5366115114
Cercon xt ML A2 disk 98 18	5366115118
Cercon xt ML A2 disk 98 25	5366115125
Cercon xt ML A3 disk 98 14	5366115214
Cercon xt ML A3 disk 98 18	5366115218
Cercon xt ML A3 disk 98 25	5366115225
Cercon xt ML A3.5 disk 98 14	5366115314
Cercon xt ML A3.5 disk 98 18	5366115318
Cercon xt ML A3.5 disk 98 25	5366115325
Cercon xt ML A4 disk 98 14	5366115414
Cercon xt ML A4 disk 98 18	5366115418
Cercon xt ML A4 disk 98 25	5366115425
Cercon xt ML B1 disk 98 14	5366115514
Cercon xt ML B1 disk 98 18	5366115518
Cercon xt ML B1 disk 98 25	5366115525
Cercon xt ML B2 disk 98 14	5366115614
Cercon xt ML B2 disk 98 18	5366115618
Cercon xt ML B2 disk 98 25	5366115625
Cercon xt ML B3 disk 98 14	5366115714
Cercon xt ML B3 disk 98 18	5366115718
Cercon xt ML B3 disk 98 25	5366115725
Cercon xt ML B4 disk 98 14	5366115814
Cercon xt ML B4 disk 98 18	5366115818
Cercon xt ML B4 disk 98 25	5366115825
Cercon xt ML C1 disk 98 14	5366115914
Cercon xt ML C1 disk 98 18	5366115918
Cercon xt ML C1 disk 98 25	5366115925
Cercon xt ML C2 disk 98 14	5366116014
Cercon xt ML C2 disk 98 18	5366116018
Cercon xt ML C2 disk 98 25	5366116025
Cercon xt ML C3 disk 98 14	5366116114
Cercon xt ML C3 disk 98 18	5366116118
Cercon xt ML C3 disk 98 25	5366116125

Cercon xt ML C4 disk 98 14	5366116214
Cercon xt ML C4 disk 98 18	5366116218
Cercon xt ML C4 disk 98 25	5366116225
Cercon xt ML D2 disk 98 14	5366116314
Cercon xt ML D2 disk 98 18	5366116318
Cercon xt ML D2 disk 98 25	5366116325
Cercon xt ML D3 disk 98 14	5366116414
Cercon xt ML D3 disk 98 18	5366116418
Cercon xt ML D3 disk 98 25	5366116425
Cercon xt ML D4 disk 98 14	5366116514
Cercon xt ML D4 disk 98 18	5366116518
Cercon xt ML D4 disk 98 25	5366116525
Cercon xt ML BL2 disk 98 14	5366116714
Cercon xt ML BL2 disk 98 18	5366116718
Cercon xt ML BL2 disk 98 25	5366116725



Cercon®ht ML

High translucent multilayer zirconia

Aesthetics meets Strength

- **Strength & Safety** – From single crowns up to 14-unit bridges, based on the high bending strength (gradient from 1200 MPa in the dentin up to 750 MPa in the incisal) and over 20 years of documented Cercon experience.
- **Natural Aesthetics** – The gradient of the layers matching natural tooth aesthetics in one disc without the need to veneer or stain.
- **Approved Color Accuracy*** – Excellent color matching for all 16 VITA** shades and BL2.
- **Easy to use** – Easy to nest and same milling and firing programs for all Cercon zirconia discs.

As a Team: Cercon ht ML and Cercon xt ML are the Zirconia Multilayer-Solution.

dentsplysirona.com

THE DENTAL
SOLUTIONS
COMPANY™



* Based on our trusted True Color Technology
** VITA is a registered brand of VITA Zahnfabrik

Cercon® ht ML– Disc 98 mm

Cercon extra translucent zirconia for enhanced aesthetics



Cercon ht ML A1 disk 98 14	5366095014
Cercon ht ML A1 disk 98 18	5366095018
Cercon ht ML A1 disk 98 25	5366095025
Cercon ht ML A2 disk 98 14	5366095114
Cercon ht ML A2 disk 98 18	5366095118
Cercon ht ML A2 disk 98 25	5366095125
Cercon ht ML A3 disk 98 14	5366095214
Cercon ht ML A3 disk 98 18	5366095218
Cercon ht ML A3 disk 98 25	5366095225
Cercon ht ML A3.5 disk 98 14	5366095314
Cercon ht ML A3.5 disk 98 18	5366095318
Cercon ht ML A3.5 disk 98 25	5366095325
Cercon ht ML A4 disk 98 14	5366095414
Cercon ht ML A4 disk 98 18	5366095418
Cercon ht ML A4 disk 98 25	5366095425
Cercon ht ML B1 disk 98 14	5366095514
Cercon ht ML B1 disk 98 18	5366095518
Cercon ht ML B1 disk 98 25	5366095525
Cercon ht ML B2 disk 98 14	5366095614
Cercon ht ML B2 disk 98 18	5366095618
Cercon ht ML B2 disk 98 25	5366095625
Cercon ht ML B3 disk 98 14	5366095714
Cercon ht ML B3 disk 98 18	5366095718
Cercon ht ML B3 disk 98 25	5366095725
Cercon ht ML B4 disk 98 14	5366095814
Cercon ht ML B4 disk 98 18	5366095818
Cercon ht ML B4 disk 98 25	5366095825
Cercon ht ML C1 disk 98 14	5366095914
Cercon ht ML C1 disk 98 18	5366095918
Cercon ht ML C1 disk 98 25	5366095925
Cercon ht ML C2 disk 98 14	5366096014
Cercon ht ML C2 disk 98 18	5366096018
Cercon ht ML C2 disk 98 25	5366096025
Cercon ht ML C3 disk 98 14	5366096114
Cercon ht ML C3 disk 98 18	5366096118
Cercon ht ML C3 disk 98 25	5366096125

Cercon ht ML C4 disk 98 14	5366096214
Cercon ht ML C4 disk 98 18	5366096218
Cercon ht ML C4 disk 98 25	5366096225
Cercon ht ML D2 disk 98 14	5366096314
Cercon ht ML D2 disk 98 18	5366096318
Cercon ht ML D2 disk 98 25	5366096325
Cercon ht ML D3 disk 98 14	5366096414
Cercon ht ML D3 disk 98 18	5366096418
Cercon ht ML D3 disk 98 25	5366096425
Cercon ht ML D4 disk 98 14	5366096514
Cercon ht ML D4 disk 98 18	5366096518
Cercon ht ML D4 disk 98 25	5366096525
Cercon ht ML BL2 disk 98 14	5366096714
Cercon ht ML BL2 disk 98 18	5366096718
Cercon ht ML BL2 disk 98 25	5366096725

Lucitone

Denture base resin



Dg-Lucitone Liq.(430 MI)	684309
Dg-Lucitone Liq.(4 Bt X 946MI)	684315
Dg-Lucitone 199 Std Org. 120MI	688103
Dg-Lucitone 199 Orig.(3X430MI)	688105
Lucitone 199 Pdr Org. 630G/Btl	688106
Lucitone 199 25Lb Pwdr (Orig.)	688111
Lucitone 199 108-U Pdr Org/5Lb	688120
Lucitone 199 630 Gm Light	688206
Lucitone 199 Pdr Lt Rd.Pnk630G	688306
Lucitone 199 Pdr Lt Rd.Pnk25Lb	688311
Lucitone 199 Pdr Dk Pnk/ 630G	688406



Celtra® Press System

Developed to make a difference

Ordering Information



Celtra® Kits

Celtra Press Starter Kit	5365400113
Celtra Press Starter Kit BL2/A2	5365400114
Celtra Ceram Basic Kit	600601
Celtra Ceram Trial Kit	600602
Celtra Ceram Esthetic Kit	600603
Celtra Ceram Gingiva Kit	600604
Celtra Ceram Shade Series Kit	600605
Dentsply Sirona Universal Stain Kit	600700
Celtra TwisTec Finishing Set	9 PCS 5359909112

Celtra® Press Ingots

MT/LT BL1	3x6 G	5365400166
MT/LT BL2	3x6 G	5365400176
LT A1	3x6 G	5365400006
LT A2	3x6 G	5365400016
LT A3	3x6 G	5365400026
LT B1	3x6 G	5365400056
LT B3	3x6 G	5365400076
LT C1	3x6 G	5365400096
LT C3	3x6 G	5365400116
LT D2	3x6 G	5365400136
LT D3	3x6 G	5365400146
MT A1	3x6 G	5365400206
MT A2	3x6 G	5365400216
MT A3	3x6 G	5365400226
MT B1	3x6 G	5365400256
MT B3	3x6 G	5365400276
MT C1	3x6 G	5365400296
MT C3	3x6 G	5365400316
MT D2	3x6 G	5365400336
MT D3	3x6 G	5365400346
HT i1	3x6 G	5365400526
HT i2	3x6 G	5365400536
HT i3	3x6 G	5365400546
MO 1	3x6 G	5365400416
MO 2	3x6 G	5365400426
MT/LT BL1	5x3 G	5365400167
MT/LT BL2	5x3 G	5365400177
LT A1	5x3 G	5365400007
LT A2	5x3 G	5365400017
LT A3	5x3 G	5365400027
LT B1	5x3 G	5365400057
LT B3	5x3 G	5365400077
LT C1	5x3 G	5365400097
LT C3	5x3 G	5365400117
LT D2	5x3 G	5365400137
LT D3	5x3 G	5365400147

MT A1	5x3 G	5365400207
MT A2	5x3 G	5365400217
MT A3	5x3 G	5365400227
MT B1	5x3 G	5365400257
MT B3	5x3 G	5365400277
MT C1	5x3 G	5365400297
MT C3	5x3 G	5365400317
MT D2	5x3 G	5365400337
MT D3	5x3 G	5365400347
HT i1	5x3 G	5365400527
HT i2	5x3 G	5365400537
HT i3	5x3 G	5365400547
MO 1	5x3 G	5365400417
MO 2	5x3 G	5365400427

Investment

Celtra Press Investment	45x100 G	5365400166
Celtra Press Investment Liquid	1000 ML	5325400320
Ring Former for 100 g Muffle		5365901206
Ring Former for 200 g Muffle		5365901207
Sprue Former for 100 g Muffle		5365901208
Sprue Former for 200 g Muffle		5365901209
Celtra Press Investment Plunger	25 PCS.	5365901211



Accessories

Celtra Undercoating wax	60 G	5 3 2 3 5 1 0 0 1 1
Celtra Modeling wax	60 G	5 3 2 3 5 1 0 0 1 9
CC Shade Guide - Basic Indicator 1		601002
CC Shade Guide - Basic Indicator 2		601003
CC Shade Guide - Esthetic Indicator		601004
CC Shade Guide - Gingiva Indicator		601005
DS Shade Guide - Stains Indicator		601006
DS Modeling Liquid DU	15 ML	615126
DS Modeling Liquid DU	50 ML	650126
DS Stain & Glaze Liquid	15 ML	601315
DS Stain & Glaze Liqui	50 ML	601350
Firing Pad	3 PCS	5365901205

Celtra® Ceram porcelain

Dentin 101 - 1M1*	15 G	615700
Dentin 102 - 1M2*	15 G	615701
Dentin 203 - 2L1.5*	15 G	615702
Dentin 204 - 2L2.5*	15 G	615703
Dentin 205 - 2M1*	15 G	615704
Dentin 206 - 2M2*	15 G	615705
Dentin 207 - 2M3*	15 G	615706
Dentin 208 - 2R1.5*	15 G	615707
Dentin 209 - 2R2.5*	15 G	615708
Dentin 310 - 3L1.5*	15 G	615709
Dentin 311 - 3L2.5*	15 G	615710
Dentin 312 - 3M1*	15 G	615711
Dentin 313 - 3M2*	15 G	615712
Dentin 314 - 3M3*	15 G	615713
Dentin 315 - 3R1.5*	15 G	615714
Dentin 316 - 3R2.5*	15 G	615715
Dentin 417 - 4L1.5*	15 G	615716
Dentin 418 - 4L2.5*	15 G	615717
Dentin 419 - 4M1*	15 G	615718
Dentin 420 - 4M2*	15 G	615719
Dentin 421 - 4M3*	15 G	615720
Dentin 422 - 4R1.5*	15 G	615721
Dentin 423 - 4R2.5*	15 G	615722
Dentin 524 - 5M1*	15 G	615723
Dentin 525 - 5M2*	15 G	615724
Dentin 526 - 5M3*	15 G	615725
Dentin - BL1	15 G	615130
Dentin - BL2	15 G	615131
Dentin - BL3	15 G	615132
Dentin - BL4	15 G	615133

Dentin - A1	15 G	615134
Dentin - A2	15 G	615135
Dentin - A3	15 G	615136
Dentin - A3.5	15 G	615137
Dentin - A4	15 G	615138
Dentin - B1	15 G	615139
Dentin - B2	15 G	615140
Dentin - B3	15 G	615141
Dentin - B4	15 G	615142
Dentin - C1	15 G	615143
Dentin - C2	15 G	615144
Dentin - C3	15 G	615145
Dentin - C4	15 G	615146
Dentin - D2	15 G	615147
Dentin - D3	15 G	615148
Dentin - D4	15 G	615149
Opaceous Dentin - OD0	15 G	615150
Opaceous Dentin - OD1	15 G	615151
Opaceous Dentin - OD2	15 G	615152
Opaceous Dentin - OD3	15 G	615153
Opaceous Dentin - OD4	15 G	615154
Opaceous Dentin - OD5	15 G	615155
Opaceous Dentin - OD6	15 G	615156
Dentin Effect - DE1, Blue	15 G	615161
Dentin Effect - DE2, Green-brown	15 G	615162
Dentin Effect - DE3, Grey	15 G	615163
Dentin Effect - DE4, Orange-brown	15 G	615164
Dentin Effect - DE5, Pink	15 G	615165
Dentin Effect - DE6, Violet	15 G	615166
Dentin Effect - DE7, Yellow	15 G	615167
Dentin Effect - DE8, White	15 G	615168
Dentin Effect - DE9, Orange	15 G	615169
Dentin Gingiva - DG1, Pink	15 G	615171
Dentin Gingiva - DG2, Reddish-pink	15 G	615172
Dentin Gingiva - DG3, Salmon	15 G	615173
Dentin Gingiva - DG4, Dark	15 G	615174
Dentin Gingiva - DG5, Mahogany	15 G	615175
Power Dentin - PD1	15 G	615181
Power Dentin - PD2	15 G	615182
Power Dentin - PD3	15 G	615183
Power Dentin - PD4	15 G	615184
Power Dentin - PD5	15 G	615185
Power Dentin - PD6	15 G	615186
Enamel - E1, Extra-light	15 G	615201
Enamel - E2, Light	15 G	615102
Enamel - E3, Medium	15 G	615103
Enamel - E4, Transparent	15 G	615204
Enamel - E5, White	15 G	615205

Enamel - E6, Dark	15 G	615206
Enamel Opal - EO1, Extra-light	15 G	615211
Enamel Opal - EO2, Light	15 G	615212
Enamel Opal - EO3, Medium	15 G	615213
Enamel Opal - EO4, Transparent	15 G	615214
Enamel Opal - EO5, HT	15 G	615215
Enamel Opal - EO6, LT	15 G	615216
Enamel Effect - EE1, Sunrise	15 G	615221
Enamel Effect - EE2, Violet	15 G	615222
Enamel Effect - EE3, Sunset	15 G	615223
Enamel Effect - EE4, Fog	15 G	615224
Enamel Effect - EE5, Sky	15 G	615225
Enamel Effect - EE6, Ivory	15 G	615226
Add-on Correction - C1, Light	15 G	615401
Add-on Correction - C2, Medium	15 G	615402
Add-on Correction - C3, Dark	15 G	615403
Add-on Correction - C4, Transparent	15 G	615404
Add-on Gingiva - G1, Pink	15 G	615411
Add-on Gingiva - G2, Reddish-pink	15 G	615412
Add-on Gingiva - G3, Salmo	15 G	615413
Add-on Gingiva - G4, Dark	15 G	615414
Add-on Gingiva - G5, Mahogany	15 G	615415

Stains and Overglaz

DS Universal Body Stain - S0	5 G	605520
DS Universal Body Stain - S1	5 G	605521
DS Universal Body Stain - S2	5 G	605522
DS Universal Body Stain - S3	5 G	605523
DS Universal Body Stain - S4	5 G	605524
DS Universal Incisal Stain - i1	5 G	605531
DS Universal Incisal Stain - i2	5 G	605532
DS Universal Stain - White	5 G	605500
DS Universal Stain - Crème	5 G	605501
DS Universal Stain - Sunset	5 G	605502
DS Universal Stain - Copper	5 G	605503
DS Universal Stain - Khaki	5 G	605504
DS Universal Stain - Olive	5 G	605505
DS Universal Stain - Mahoga	5 G	605506
DS Universal Stain - Violet	5 G	605507
DS Universal Stain - Raspberry	5 G	605508
DS Universal Stain - Purple	5 G	605509
DS Universal Stain - Grey	5 G	605510
DS Universal Stain - Chestnut	5 G	605511
DS Universal Stain - Blue	5 G	605512
DS Universal Stain - Pink	5 G	605513
DS Universal Overglaze	5 G	605540
DS Universal Overglaze - High Flu	5 G	605542




Overview Celtra® Kits



Convenient kit that includes everything necessary to enjoy the benefits of this smart, integrated product system. Includes ingots, investment and ceramics as well as several key accessories to maximize performance.



Offers convenience of most popular shades, ancillary porcelains and accessories to quickly get started enjoying the benefits of this exceptional veneering system.



Mini-kit offers convenience for the esthetic-minded technician wishing to create individualized characterization of Celtra Press and lithium disilicate frameworks.

Celtra® Press Starter Kit REF 5365400113

- **Celtra Press** | 5x3G
MT A2, MT A3
- **DENTSPLY Prosthetics Die Material** | 4G
F10, F12
- Celtra Undercoating Wax | 60G
- Celtra Modeling Wax | 60G
- Celtra Press Investment | 2 x 100G
- Celtra Press Investment Liquid | 100ML
- Celtra Press Crucible former | 100G
- Celtra Press Ring Former | 100G
- Celtra Press Investment Plunger | 2 PCS
- Firing Pad | 3 PCS
- **Celtra Ceram Enamel** | 15G
E1, Extra-light
E2, Light
- DS Universal Overglaze, High Flu | 5G
- Universal Stain & Glaze Liquid | 15ML
- DS Modeling Liquid DU | 15ML

Celtra® Ceram Basic Kit REF 600601

- **Celtra Ceram Dentin** | 15G
BL2, A1, A2, A3, A3.5, B1, B2, B3, C1, C2, C3, D2, D3
- **Celtra Ceram Enamel** | 15G
E1, Extra-light
E2, Light
E3, Medium
E4, Transparent
E5, White
- **Celtra Ceram Enamel Opal** | 15G
EO1, Extra-light
EO2, Light
EO3, Medium
EO4, Transparent
EO5, HT
EO6, LT
- **Celtra Ceram Add-on Correction** | 15G
C1, Light
C2, Medium
C3, Dark
C4, Transparent
- DS Universal Overglaze, High Flu | 5G
- DS Modeling Liquid DU | 50ML
- Universal Stain & Glaze Liquid | 50ML
- DS Portioner
- MagicWizard Brush
- MagicStain Brush
- OccluSpeed Brush
- Shade Guide - Basic Indicator 1
- Shade Guide - Basic Indicator 2

Celtra® Ceram Aesthetic Kit REF 600603

- **Celtra Ceram Dentin Effect** | 15G
DE1, Blue
DE2, Green-brown
DE3, Grey
DE4, Orange-brown
DE5, Pink
DE6, Violet
DE7, Yellow
DE8, White
DE9, Orange
- **Celtra Ceram Power Dentin** | 15G
PD1, PD2, PD3, PD4, PD5, PD6
- **Celtra Ceram Opaceous Dentin** | 15G
OD0 - OD6
- **Celtra Ceram Enamel Effect** | 15G
EE1, Sunrise
EE2, Violet
EE3, Sunset
EE4, Fog
EE5, Sky
EE6, Ivory
- **Celtra Ceram Dentin Gingiva** | 15G
DG3, Salmon
- **Celtra Ceram Add-on Gingiva** | 15G
G3, Salmon
- Shade Guide - Esthetic Indicator
- Shade Guide - Gingiva Indicator



Convenient kit promoting stain and overglaze necessary for achieving 16 VITA Classical shades for Cercon, Celtra Press and lithium disilicate frameworks. Wide compatibility and exceptional ease of use make these items a must have in any dental laboratory seeking efficiency.

Dentsply Sirona Universal Stain Kit REF 600700

- **DS Universal Body Stain** | 5G
S0, S1, S2, S3, S4
- **DS Universal Incisal Stain** | 5G
i1, i2
- **DS Universal Stain** | 5G
White, Creme, Sunset, Copper, Khaki, Olive, Mahogany, Violet, Purple, Raspberry, Grey, Chestnut, Blue, Pink
- Universal Stain & Glaze Liquid | 15ML
- DS Universal Overglaze | 5G
- DS Universal Overglaze, High Flu | 5G
- Shade Guide - Stains Indicator
- MagicStain Brush



Mini-kit with the necessary components to create life-like gingiva effects, especially for implant-retained prostheses.

Celtra® Ceram Gingiva Kit REF 600604

- **Celtra Ceram Dentin Gingiva** | 15G
DG1, Pink
DG2, Reddish Pink
DG3, Salmon
DG4, Dark
DG5, Mahogany
- **Celtra Ceram Add-on Gingiva** | 15G
G1, Pink
G2, Reddish Pink
G3, Salmon
G4, Dark
G5, Mahogany
- DS Modeling Liquid DU | 50ML
- Shade Guide - Gingiva Indicator



Just the basics for experiencing the incredible consistency, robustness and esthetic capabilities of the Celtra Ceram system.

Celtra® Ceram Trial Kit REF 600602

- **Celtra Ceram Dentin** | 15G
A2, A3
- **Celtra Ceram Enamel** | 15G
E1, Extra-light
E2, Light
E4, Transparent
- **Celtra Ceram Add-on Correction** | 15G
C3, Dark
- DS Universal Overglaze, High Flu | 5G
- DS Modeling Liquid DU | 15 ML
- Universal Stain & Glaze Liquid | 15ML



Celtra® Ceram Shade Series Kit REF 600605

This Kit offers the basics for experiencing the incredible 3D* components and a darker enamel shade.

- **Celtra Ceram Dentin** | 15G
Dentin 101, Dentin 102, Dentin 203, Dentin 204, Dentin 205, Dentin 206, Dentin 207, Dentin 208, Dentin 209, Dentin 310, Dentin 311, Dentin 312, Dentin 313, Dentin 314, Dentin 315, Dentin 316, Dentin 417, Dentin 418, Dentin 419, Dentin 420, Dentin 421, Dentin 422, Dentin 423, Dentin 524, Dentin 525, Dentin 526
- **Celtra Ceram Enamel** | 15G
E6, Dark

* VITA, VITA classical A1-D4, and 3D-MASTER are registered trademarks of VITA Zahnfabrik H. Rauter GmbH & Co. KG. The product shades referred to in our comparative list are Dentsply Sirona's closest approximation to the referenced VITA 3D-Master shades.



Celtra® Ceram
All-Ceramic Veneering Porcelain

Versatility that delivers

THE DENTAL
SOLUTIONS
COMPANY™

 **Dentsply
Sirona**

Celtra® Ceram

Versatility that delivers

Celtra Ceram introduces the versatility of a single veneering system designed exclusively for all-ceramic substructures including lithium silicates and zirconia oxides.

- + Versatile concept
- + Premium aesthetics
- + Consistent performance

Celtra Ceram is suitable for veneering all-ceramic frameworks.
Coefficient of Thermal Expansion (CTE): $9.0 \times 10^{-6} \text{ K}^{-1}$ (25-500°C)

- + Celtra Press, Zirconia-reinforced Lithium Silicate (ZLS) frameworks: $9.7 \times 10^{-6} \text{ K}^{-1}$ (CTE 25-500°C)
- + Lithium Disilicate: $10.0 - 10.5 \times 10^{-6} \text{ K}^{-1}$ (25-500°C)
- + Cercon® ht frameworks / Cercon Base: $10.5 \times 10^{-6} \text{ K}^{-1}$ (25-500°C)
- + Cercon xt frameworks: $10.1 \times 10^{-6} \text{ K}^{-1}$ (25-500°C)
- + Zirconia frameworks: $10.1 - 11.0 \times 10^{-6} \text{ K}^{-1}$ (25-500°C)

EXCEPTIONAL PERFORMANCE

For incredibly stable results



A proprietary feldspathic base formulation delivers an incredibly stable and robust result – every time, even after multiple firings

SHADE CONFIDENCE

For lifelike, reproducible veneerings



Out-of-the-bottle shade accuracy: Fewer remakes, improved consistency





Enamel

Enamel - E1, Extra-light	15g	615201
Enamel - E2, Light	15g	615202
Enamel - E3, Medium	15g	615203
Enamel - E4, Transparent	15g	615204
Enamel - E5, White	15g	615205
Enamel - E6, Dark	15g	615206

Enamel Opal

Enamel Opal - EO1, Extra-light	15g	615211
Enamel Opal - EO2, Light	15g	615212
Enamel Opal - EO3, Medium	15g	615213
Enamel Opal - EO4, Transparent	15g	615214
Enamel Opal - EO5, HT	15g	615215
Enamel Opal - EO6, LT	15g	615216

Enamel Effect

Enamel Effect - EE1, Sunrise	15g	615221
Enamel Effect - EE2, Violet	15g	615222
Enamel Effect - EE3, Sunse	15g	615223
Enamel Effect - EE4, Fog	15g	615224
Enamel Effect - EE5, Sky	15g	615225
Enamel Effect - EE6, Ivory	15g	615226

Add-On Correction

Add-on Correction - C1, Light	15g	615401
Add-on Correction - C2, Medium	15g	615402
Add-on Correction - C3, Dark	15g	615403
Add-on Correction - C4, Transparent	15g	615404

Add-On Gingiva

Add-on Gingiva - G1, Pink	15g	615411
Add-on Gingiva - G2, Reddish-pink	15g	615412
Add-on Gingiva - G3, Salmon	15g	615413
Add-on Gingiva - G4, Dark	15g	615414
Add-on Gingiva - G5, Mahogany	15g	615415

Kit Configurations

Basic Kit (BL2, A1-A3.5, B1-B3, C1-C3, D2-D3)	EA	600601
CC Kit - Trial (A2, A3)	EA	600602
CC Kit - Aesthetic	EA	600603
CC Kit - Gingiva	EA	600604
Shade Series Kit	EA	600605

Accessories

Shade Indicator - Basic 1	EA	D601002
Shade Indicator - Basic 2	EA	D601003
Shade Indicator - Aesthetic	EA	D601004
Shade Indicator - Gingiva	EA	D601005
Dentsply Sirona Modeling Liquid U	100mL	301126
Dentsply Sirona Modeling Liquid E	100mL	301127
Dentsply Sirona Modeling Liquid DU	50mL	650126
Ducera Liquid SD	50mL	5368272054
Ducera Liquid FORM	50mL	5368272034
Dentsply Sirona Stain and Glaze Liquid	15mL	601315

Dentsply Sirona Universal Stain & Overglaze

DS Universal Stain Kit (all items below)	EA	600700
DS Universal Body Stain - S0	5g	605520
DS Universal Body Stain - S1	5g	605521
DS Universal Body Stain - S2	5g	605522
DS Universal Body Stain - S3	5g	605523
DS Universal Body Stain - S4	5g	605524
DS Universal Incisal Stain - i1	5g	605531
DS Universal Incisal Stain - i2	5g	605532
DS Universal Stain - White	5g	605500
DS Universal Stain - Crème	5g	605501
DS Universal Stain - Sunset	5g	605502
DS Universal Stain - Copper	5g	605503
DS Universal Stain - Khaki	5g	605504
DS Universal Stain - Olive	5g	605505
DS Universal Stain - Mahogany	5g	605506
DS Universal Stain - Violet	5g	605507
DS Universal Stain - Raspberry	5g	605508
DS Universal Stain - Purple	5g	605509
DS Universal Stain - Grey	5g	605510
DS Universal Stain - Chestnut	5g	605511
DS Universal Stain - Blue	5g	605512
DS Universal Stain - Pink	5g	605513
DS Universal Overglaze	5g	605540
DS Universal Overglaze - High Flu	5g	605541

Dentin

Dentin - BL1	15g	615130
Dentin - BL2	15g	615131
Dentin - BL3	15g	615132
Dentin - BL4	15g	615133
Dentin - A1	15g	615134
Dentin - A2	15g	615135
Dentin - A3	15g	615136
Dentin - A3.5	15g	615137
Dentin - A4	15g	615138
Dentin - B1	15g	615139
Dentin - B2	15g	615140
Dentin - B3	15g	615141
Dentin - B4	15g	615142
Dentin - C1	15g	615143
Dentin - C2	15g	615144
Dentin - C3	15g	615145
Dentin - C4	15g	615146
Dentin - D2	15g	615147
Dentin - D3	15g	615148
Dentin - D4	15g	615149

Opaceous Dentin

Opaceous Dentin - OD0	15g	615150
Opaceous Dentin - OD1	15g	615151
Opaceous Dentin - OD2	15g	615152
Opaceous Dentin - OD3	15g	615153
Opaceous Dentin - OD4	15g	615154
Opaceous Dentin - OD5	15g	615155
Opaceous Dentin - OD6	15g	615156

Dentin Effect

Dentin Effect - DE1, Blue	15g	615161
Dentin Effect - DE2, Green-brown	15g	615162
Dentin Effect - DE3, Grey	15g	615163
Dentin Effect - DE4, Orange-brown	15g	615164
Dentin Effect - DE5, Pink	15g	615165
Dentin Effect - DE6, Violet	15g	615166
Dentin Effect - DE7, Yellow	15g	615167
Dentin Effect - DE8, White	15g	615168
Dentin Effect - DE9, Orange	15g	615169

Dentin Gingiva

Dentin Gingiva - DG1, Pink	15g	615171
Dentin Gingiva - DG2, Reddish-pink	15g	615172
Dentin Gingiva - DG3, Salmon	15g	615173
Dentin Gingiva - DG4, Dark	15g	615174
Dentin Gingiva - DG5, Mahogany	15g	615175

Power Dentin

Power Dentin - PD1	15g	615181
Power Dentin - PD2	15g	615182
Power Dentin - PD3	15g	615183
Power Dentin - PD4	15g	615184
Power Dentin - PD5	15g	615185
Power Dentin - PD6	15g	615186

Celtra Ceram, Shade Series Kit <i>Includes all items below</i>			600605
Dentin 101	1M1	15g	615700
Dentin 102	1M2	15g	615701
Dentin 203	2L1.5	15g	615702
Dentin 204	2L2.5	15g	615703
Dentin 205	2M1	15g	615704
Dentin 206	2M2	15g	615705
Dentin 207	2M3	15g	615706
Dentin 208	2R1.5	15g	615707
Dentin 209	2R2.5	15g	615708
Dentin 310	3L1.5	15g	615709
Dentin 311	3L2.5	15g	615710
Dentin 312	3M1	15g	615711
Dentin 313	3M2	15g	615712
Dentin 314	3M3	15g	615713
Dentin 315	3R1.5	15g	615714
Dentin 316	3R2.5	15g	615715
Dentin 417	4L1.5	15g	615716
Dentin 418	4L2.5	15g	615717
Dentin 419	4M1	15g	615718
Dentin 420	4M2	15g	615719
Dentin 421	4M3	15g	615720
Dentin 422	4R1.5	15g	615721
Dentin 423	4R2.5	15g	615722
Dentin 524	5M1	15g	615723
Dentin 525	5M2	15g	615724
Dentin 526	5M3	15g	615725



Ceramco® 3 veneering ceramics

Ceramco 3 is our value-driven, high-fusing veneering system that introduces a simple yet highly efficient layering concept accommodating a wide array of technique and proficiency levels in the dental laboratory. Its patented color concept and proprietary mineral base allow for better individual control, so technicians have the freedom to work the way they prefer to work and always achieve predictable results.



Photo courtesy Angelo Gallochio, Dental technician / Olivieri Studio Dentistico, Borgomanero, Italy

Simple concept

- + Complete system providing stability and freedom so technicians can work how they prefer and achieve a great outcome

Easy to use

- + Exceptional ease of use for all skill levels using simple, color-coded system satisfying basic and advanced techniques

Consistent performance

- + Shade accurate, stable and forgiveness achieved using a proprietary mineral source and a patented¹ color control concept
- + Recently enhanced Dentin Modifiers deliver chromatic and hue effects to mimic natural variation

¹US Patent No. 5,482,732

Consistent, aesthetic results from every technique in your lab, every time

With over 166 million crowns made in laboratories worldwide, Ceramco® 3 is a proven veneering system that delivers the results you demand every time.

- + Simple concept
- + Easy to use
- + Consistent performance

The Ceramco 3 veneering porcelains may be used for all single and multiple unit porcelain fused-to-metal fixed prosthodontics.



Ceramco 3 is great because of its reliability, the translucency of the dentine and the wide variety of enamels Ceramco 3 has. It guarantees a perfect camouflage in every situation.



Pierpaolo Mola
Odontomaster
Italy

166,329,226 Ceramco® 3 crowns created by laboratories worldwide¹



¹and counting, using Ceramco 3 since the product introduction.
Photos courtesy of Pierpaolo Mola, Odontomaster, Italy



Powder Opaqu 2mL

Paste Opaque - Base Paste	301040
Paste Opaque - Base Paste Illuminé	301092
Paste Opaque - Illuminé i1/i3	301044
Paste Opaque - Illuminé i2/4	301053
Paste Opaque - Illuminé i5/i7	301042
Paste Opaque - Illuminé i6/i8	301051
Paste Opaque - A1	301045
Paste Opaque - A2	301046
Paste Opaque - A3	301047
Paste Opaque - A3.5	301048
Paste Opaque - A4	301049
Paste Opaque - B1	301054
Paste Opaque - B2	301055
Paste Opaque - B3	301056
Paste Opaque - B4	301057
Paste Opaque - C1	301058
Paste Opaque - C2	301059
Paste Opaque - C3	301060
Paste Opaque - C4	301061
Paste Opaque - D2	301062
Paste Opaque - D3	301064
Paste Opaque - D4	301065
Paste Opaque - Shade Series 101	303010
Paste Opaque - Shade Series 102	303011
Paste Opaque - Shade Series 203	303012
Paste Opaque - Shade Series 204	303013
Paste Opaque - Shade Series 205	303014
Paste Opaque - Shade Series 206	303015
Paste Opaque - Shade Series 207	303016
Paste Opaque - Shade Series 208	303017
Paste Opaque - Shade Series 209	303018
Paste Opaque - Shade Series 310	303019
Paste Opaque - Shade Series 311	303020
Paste Opaque - Shade Series 312	303021
Paste Opaque - Shade Series 313	303022
Paste Opaque - Shade Series 314	303023
Paste Opaque - Shade Series 315	303024
Paste Opaque - Shade Series 316	303025
Paste Opaque - Shade Series 417	303026
Paste Opaque - Shade Series 418	303027
Paste Opaque - Shade Series 419	303028
Paste Opaque - Shade Series 420	303029
Paste Opaque - Shade Series 421	303030
Paste Opaque - Shade Series 422	303031
Paste Opaque - Shade Series 423	303032
Paste Opaque - Shade Series 524	303033
Paste Opaque - Shade Series 525	303034
Paste Opaque - Shade Series 526	303035

Paste Opaque Crystals 10g

Opaque Crystals, White	301080
Opaque Crystals, Light	301081
Opaque Crystals, Medium	301082
Opaque Crystals, Dark	301083

Opaque Modifier 3.5mL / 28g

	3.5mL	28g
Opaque Modifier, White	301070	301030
Opaque Modifier, Yellow	301071	301031
Opaque Modifier, Sienna	301072	301032
Opaque Modifier, Blue	301073	
Opaque Modifier, Violet		301033
Opaque Modifier, Gray	301074	301034
Opaque Modifier, Tan	301075	301035
Opaque Modifier, Pink	301076	301036
Opaque Modifier, Ochre	301077	301037
Opaque Modifier, Orange	301078	301038

Opaque Corrector

	29.4g
Opaque Corrector, White	301413
Opaque Corrector, Light	301410
Opaque Corrector, Medium	301411
Opaque Corrector, Dark	301412

Dentin

	10g	28g	113g
Dentin - Illuminé i1	305133	301133	304130
Dentin - Illuminé i2	305142	301142	304131
Dentin - Illuminé i3	305132	301132	304132
Dentin - Illuminé i4	305141	301141	304133
Dentin - Illuminé i5	305131	301131	304134
Dentin - Illuminé i6	305140	301140	304135
Dentin - Illuminé i7	305130	301130	304136
Dentin - Illuminé i8	305139	301139	304137
Dentin - A1	305134	301134	304138
Dentin - A2	305135	301135	304139
Dentin - A3	305136	301136	304140
Dentin - A3.5	305137	301137	304141
Dentin - A4	305138	301138	304142
Dentin - B1	305143	301143	304143
Dentin - B2	305144	301144	304144
Dentin - B3	305145	301145	304145

	10g	28g	113g
Dentin - B4	305146	301146	304146
Dentin - C1	305147	301147	304147
Dentin - C2	305148	301148	304148
Dentin - C3	305149	301149	304149
Dentin - C4	305150	301150	304150
Dentin - D2	305151	301151	304151
Dentin - D3	305152	301152	304152
Dentin - D4	305153	301153	304153
Dentin - 101	303700	303701	
Dentin - 102	303710	303711	
Dentin - 203	303720	303721	
Dentin - 204	303730	303731	
Dentin - 205	303740	303741	
Dentin - 206	303750	303751	
Dentin - 207	303760	303761	
Dentin - 208	303770	303771	
Dentin - 209	303780	303781	
Dentin - 310	303790	303791	
Dentin - 311	303800	303801	
Dentin - 312	303810	303811	
Dentin - 313	303820	303821	
Dentin - 314	303830	303831	
Dentin - 315	303840	303841	
Dentin - 316	303850	303851	
Dentin - 417	303860	303861	
Dentin - 418	303870	303871	
Dentin - 419	303880	303881	
Dentin - 420	303890	303891	
Dentin - 421	303900	303901	
Dentin - 422	303910	303911	
Dentin - 423	303920	303921	
Dentin - 524	303930	303931	
Dentin - 525	303940	303941	
Dentin - 526	303950	303951	



Opaceous Dentin

	10g	28g	113g
Opaceous Dentin - Illuminé i1	305103	301103	304100
Opaceous Dentin - Illuminé i2	305112	301112	304101
Opaceous Dentin - Illuminé i3	305102	301102	304102
Opaceous Dentin - Illuminé i4	305111	301111	304103
Opaceous Dentin - Illuminé i5	305101	301101	304104
Opaceous Dentin - Illuminé i6	305110	301110	304105
Opaceous Dentin - Illuminé i7	305100	301100	304106
Opaceous Dentin - Illuminé i8	305109	301109	304107
Opaceous Dentin - A1	305104	301104	304108
Opaceous Dentin - A2	305105	301105	304109
Opaceous Dentin - A3	305106	301106	304110
Opaceous Dentin - A3.5	305107	301107	304111
Opaceous Dentin - A4	305108	301108	304112
Opaceous Dentin - B1	305113	301113	304113
Opaceous Dentin - B2	305114	301114	304114
Opaceous Dentin - B3	305115	301115	304115
Opaceous Dentin - B4	305116	301116	304116
Opaceous Dentin - C1	305117	301117	304117
Opaceous Dentin - C2	305118	301118	304118
Opaceous Dentin - C3	305119	301119	304119
Opaceous Dentin - C4	305120	301120	304120
Opaceous Dentin - D2	305121	301121	304121
Opaceous Dentin - D3	305122	301122	304122
Opaceous Dentin - D4	305123	301124	304123
Opaceous Dentin - 101	303400	303401	
Opaceous Dentin - 102	303410	303411	
Opaceous Dentin - 203	303420	303421	
Opaceous Dentin - 204	303430	303431	
Opaceous Dentin - 205	303440	303441	
Opaceous Dentin - 206	303450	303451	
Opaceous Dentin - 207	303460	303461	
Opaceous Dentin - 208	303470	303471	
Opaceous Dentin - 209	303480	303481	
Opaceous Dentin - 310	303490	303491	
Opaceous Dentin - 311	303500	303501	
Opaceous Dentin - 312	303510	303511	
Opaceous Dentin - 313	303520	303521	
Opaceous Dentin - 314	303530	303531	
Opaceous Dentin - 315	303540	303541	
Opaceous Dentin - 316	303550	303551	
Opaceous Dentin - 417	303560	303561	
Opaceous Dentin - 418	303570	303571	
Opaceous Dentin - 419	303580	303581	
Opaceous Dentin - 420	303590	303591	
Opaceous Dentin - 421	303600	303601	
Opaceous Dentin - 422	303610	303611	
Opaceous Dentin - 423	303620	303621	



	10g	28g
Opaceous Dentin - 524	303630	303631
Opaceous Dentin - 525	303640	303641
Opaceous Dentin - 526	303650	303651

Dentin Modifier 10g / 28g

	10g	28g
Dentin Modifier, High Chroma HC1	305960	301960
Dentin Modifier, High Chroma HC2	305961	301961
Dentin Modifier, High Chroma HC3	305962	301962
Dentin Modifier, High Chroma HC4	305963	301963
Dentin Modifier, High Chroma HC5	305964	301964
Dentin Modifier, High Chroma HC6	305965	301965
Dentin Modifier, White	305166	301166
Dentin Modifier, Blue	305167	301167
Dentin Modifier, Yellow	305168	301168
Dentin Modifier, Violet	305169	301169
Dentin Modifier, Gray	305170	301170
Dentin Modifier, Pink	305171	301171
Dentin Modifier, Orange	305173	301173
Dentin Modifier, Neutral	305174	
Dentin Modifier, Intense Orange	305175	301175

Natural Enamel 10g / 28g / 133g

	10g	28g	133g
Natural Enamel - White	305200	301200	304200
Natural Enamel - Extra Light	305201	301201	304201
Natural Enamel - Light	305202	301202	304202
Natural Enamel - Medium	305203	301203	304203
Natural Enamel - Dark	305204	301204	304204
Natural Enamel - Ivory	305205	301205	304205
Natural Enamel - Violet	305206	301206	304206
Natural Enamel - Orange	305207	301207	304207
Natural Enamel - Blue	305208	301208	304208
Natural Enamel - Pink	305209	301209	304209
Natural Enamel - Yellow	305210	301210	304210
Natural Enamel - Clear	305211	301211	304211
Natural Enamel - Super Clear	305212	301212	

Opal Enamel 10g / 28g

	10g	28g
Opal Enamel - White	305220	301220
Opal Enamel - Light	305221	301221
Opal Enamel - Medium	305222	301222
Opal Enamel - Clear	305223	301223

White Enamel

	10g
White Enamel - Transparent	305224

Shade Series Enamel

	10g	28g
Shade Series Enamel - Light	303050	303051
Shade Series Enamel - Medium	303060	303061

Mamelon

	10g	28g
Mamelon - Red-orange	305230	301230
Mamelon - Yellow-orange	305231	301231

Add-On

	10g	28g
Add-On - Translucent	305400	301400
Add-On - White	305401	301401
Add-On - Light	305402	301402
Add-On - Medium	305403	301403
Add-On - Dark	305404	301404

Add-On Tissue Tint

	10g	28g
Tissue Tint - Pink	305405	301405
Tissue Tint - Reddish-pin	305406	301406
Tissue Tint - Salmon	305407	301407
Tissue Tint - Dark	305408	301408

Margin

	10g	28g
Margin - Illuminé i1	305900	301900
Margin - Illuminé i2	305901	301901
Margin - Illuminé i5	305902	301902
Margin - A1	305903	301903
Margin - A3	305904	301904
Margin - A4	305905	301905
Margin - B2	305906	301906
Margin - C2	305907	301907
Margin - C3	305908	301908

Margin Modifier

	10g	28g
Margin Modifier - White	305909	301909
Margin Modifier - Peach	305910	301910
Margin Modifier - Honey-yellow	305911	301911
Margin Modifier - Yellow-red	305913	301913
Margin Modifier - Orange-brown	305914	301914
Margin Modifier - Green-brown	305915	301915

Final Margin

	10g
Final Margin - Illuminé i1	301920
Final Margin - Illuminé i2	301921
Final Margin - Illuminé i5	301922
Final Margin - A1	301923
Final Margin - A3	301924
Final Margin - A4	301925
Final Margin - B2	301926
Final Margin - C2	301927
Final Margin - C3	301928

Accessories

Shade Indicator - Dentin, Bleach & A-D	EA	301522
Shade Indicator - Dentin, Shade Series	EA	301528
Shade Indicator - Enamel	EA	301523
Shade Indicator - Modifier/Mamelon/ Add-On	EA	301524
Shade Indicator - Margin/Tissue Tint	EA	301525
Shade Indicator - DS Universal Stain	EA	D601006
Ceramco 3 Opaque Modifier Liquid	15mL	301084



Ceramco 3 Margin/Opaque Corrector Liquid	15mL	301916
Dentsply Sirona Modeling Liquid U	100mL	301126
Dentsply Sirona Modeling Liquid E	100mL	301127
Ducera Liquid SD	50mL	53 6827 2054
Ducera Liquid FORM	50mL	53 6827 203
Dentsply Sirona Stain and Glaze Liquid	15mL	601315
Opaque Crystals Collection Bowl - White	EA	301090
Opaque Crystals Collection Bowl - Light	EA	301087
Opaque Crystals Collection Bowl - Medium	EA	301088
Opaque Crystals Collection Bowl - Dark	EA	301089

Dentsply Sirona Universal Stain & Overglaze

	EA	5g
DS Universal Stain Kit (all items below)	600700	
DS Universal Body Stain - S0		605520
DS Universal Body Stain - S1		605521
DS Universal Body Stain - S2		605522
DS Universal Body Stain - S3		605523
DS Universal Body Stain - S4		605524
DS Universal Incisal Stain - i1		605531
DS Universal Incisal Stain - i2		605532
DS Universal Stain - White		605500
DS Universal Stain - Crème		605501
DS Universal Stain - Sunset		605502
DS Universal Stain - Copper		605503
DS Universal Stain - Khaki		605504
DS Universal Stain - Olive		605505
DS Universal Stain - Mahogany		605506
DS Universal Stain - Violet		605507
DS Universal Stain - Raspberry		605508
DS Universal Stain - Purple		605509
DS Universal Stain - Grey		605510
DS Universal Stain - Chestnut		605511
DS Universal Stain - Blue		605512
DS Universal Stain - Pink		605513
DS Universal Overglaze		605540
DS Universal Overglaze - High Flu		605542



Universal Stain and Overglaze All your needs in one

THE DENTAL
SOLUTIONS
COMPANY™



Universal Stain and Overglaze

Diverse application possibilities

The material properties of the Universal Stains make them suitable for use in custom-staining all-ceramic framework materials as well as for intermediate layers within veneering ceramic build-ups, for example between dentin firing and correction firing.

The stains should be applied moderately and in thin layers.

Table 1: Assignment of Shade Groups

Vita Shade Group	Bleach	A	B	C	D
Universal Stain Shade	Stain 0	Stain 1	Stain 2	Stain 3	Stain 4

Table 2: Use of Incisal Stains

Stain	Tooth Shade
Incisal Stain i1	A1, A2, A3, B1, B2, B3, B4, D2, D3, Bleach
Incisal Stain i2	A3.5, A4, C1, C2, C3, C4, D4



Angelo Gallochio, Dental technician | Olivieri Studio Dentistico, Borgomanero, Italy



Body Stain

DS Universal Body Stain - S0	5g	605520
DS Universal Body Stain - S1	5g	605521
DS Universal Body Stain - S2	5g	605522
DS Universal Body Stain - S3	5g	605523
DS Universal Body Stain - S4	5g	605524

Incisal Stain

DS Universal Incisal Stain - i1	5g	605531
DS Universal Incisal Stain - i2	5g	605532

Color Stain

DS Universal Stain - White	5g	605500
DS Universal Stain - Crème	5g	605501
DS Universal Stain - Sunset	5g	605502
DS Universal Stain - Copper	5g	605503
DS Universal Stain - Khaki	5g	605504
DS Universal Stain - Olive	5g	605505
DS Universal Stain - Mahogany	5g	605506
DS Universal Stain - Violet	5g	605507
DS Universal Stain - Raspberry	5g	605508
DS Universal Stain - Purple	5g	605509
DS Universal Stain - Grey	5g	605510
DS Universal Stain - Chestnut	5g	605511
DS Universal Stain - Blue	5g	605512
DS Universal Stain - Pink	5g	605513

Overglaze

DS Universal Overglaze	5g	605540
DS Universal Overglaze - High Flu	5g	605542

Kit Configuration

DS Universal Stain Kit (all items above)	EA	600700
---	----	--------

A finishing system for all your needs

The Universal Stain and Overglaze is a highly versatile, comprehensive stain system designed with a full spectrum of color options to meet every need. Broad compatibility makes having a single finishing system a true possibility in your dental laboratory.



Ultimate efficiency

- + Universal material compatibility eliminates multiple finishing systems in your lab

Easy to use

- + Preferred gel format simplifies training and makes it easier to achieve a consistent outcome

Versatile performance

- + Able to be used internally and externally to achieve highly individualized effects

Cosmo *HXL*

Natural and functional.
Setting the tone in prosthodontics.



Setting the tone in prosthodontics

To realize this ambitious goal, we have mobilized all our resources: intensive research, valuable raw materials, up-to-date production methods, and rigour quality control.

Cosmo *HXL*
Natural and functional.



The Natural Choice

Popular moulds, shades and excellent articulation make Cosmo *HXL* the natural choice. The success of this tooth-line is based on a perfect combination of three parameters:

- Natural aesthetics
- Excellent material properties
- Easy handling

Superior HXL Material

Cosmo teeth are made of HXL (= highly cross-linked) high quality PMMA resin. This ensures outstanding wear rate figures.

Easy Handling

The logical mould coding (tapering, round, square) makes the shape selection easier and more convenient. Moreover Cosmo *HXL* posteriors are easy to set-up due to their reduced cusp angle of 20°. All Cosmo *HXL* teeth are available in the most popular V-shades. This shade selection is simplified, too. The comprehensive mould variety helps to minimize the tooth stock without losing the flexibility to fulfill every prosthetic case.

Combination Table

Cosmo <i>HXL</i> Upper Anteriors	Cosmo <i>HXL</i> Lower Anteriors	Cosmo <i>HXL</i> Posteriors
T2	L1/L3	91
T4	L1/L3	91/92
T6	L3/L5	91/93
T8	L5/L7	92/93
T10	L5/L7	92/93
T12	L5/L9	92/93
T14	L7/L9	92/93
T16	L7/L9	93/94
T18	L7/L9	93/94
T20	L9/L11/L13	94/95
T22	L13	95
R2	L1	91
R4	L5	92/93
R6	L5/L7	92/93
R8	L7	92
S2	L1/L3	91
S4	L3/L5	91/93
S6	L3/L5	93
S8	L7	93
S10	L7	93
S12	L7	93
S18	L11/L13	95

Cosmo HXL

Upper Anteriors

Available in V Shades
(A1-D4 without C1 and B1)

Mould	Tapering	Width of the complete set in mm	Length of the centrals without neck in mm	Width of the centrals in mm
T2		40.2	8.9	7.3
T4		40.6	9.5	7.7
T6		42.2	9.6	8.2
T8		43.4	10.0	8.1
T10		42.8	10.0	7.8
T12		44.7	10.7	8.6
T14		45.2	10.2	8.4
T16		45.4	10.7	8.7
T18		47.6	10.9	9.0

Mould	Tapering	Width of the complete set in mm	Length of the centrals without neck in mm	Width of the centrals in mm
T20		51.2	11.3	9.6
T22		53.2	12.2	10.0
Mould	Round	Width of the complete set in mm	Length of the centrals without neck in mm	Width of the centrals in mm
R2		40.6	9.0	7.6
R4		42.6	9.7	7.9
R6		44.5	9.6	8.2
R8		45.7	10.7	8.4

Cosmo HXL

Upper Anteriors

Available in V Shades
(A1-D4 without C1 and B1)

Mould	Square	Width of the complete set in mm	Length of the centrals without neck in mm	Width of the centrals in mm
S2		40.8	8.6	7.7
S4		43.7	9.2	8.2
S6		43.3	9.8	8.3
S8		45.3	9.4	8.8
S10		45.2	10.1	8.2
S12		46.6	10.4	8.7
S18		50.2	11.5	9.5

Cosmo HXL

Lower Anteriors

Available in V Shades
(A1-D4 without C1 and B1)

Mould		Width of the complete set in mm	Length of the centrals without neck in mm	Width of the centrals in mm
L1		30.6	8.6	4.6
L3		31.5	8.7	4.6
L5		33.1	10.0	4.9
L7		35.9	10.0	5.7
L9		38.1	10.1	5.7
L11		39.5	11.0	5.7
L13				

Cosmo HXL

Posteriors

Available in V Shades
(A1-D4 without C1 and B1)

Mould

91



28.9 8.4

30.9 7.6

92



28.7 8.5

30.7 7.9

93



29.9 8.4

32.7 7.2

94



30.5 8.8

33.3 7.4

95



33.0 9.5

35.3 8.7

Half the width of
the set in mm
Width of the first molar
in mm

Sales through:

58004/0503/GVD



Perfectly combined components

To ensure a high level of processing reliability, we offer a comprehensive prosthodontic system, all components are perfectly combined with our toothlines.

Selectaplus/Selectaplus H

A cold-cure and a heat-cure resin for complete and partial dentures, repairs and relinings. Selectaplus is characterized by high fracture resistance, colour stability and plaque resistance.

StarLoy L

Laserable cobalt-chromium alloy for cast-metal frameworks. StarLoy L features a low melting range, castings of very high quality, and is extra spring-hard.

SilberStar

Investment material for metal casting with optimized flow characteristics and a stable model surface.

ReproStar

Addition-curing two-component silicone, suitable for all duplicating techniques. Stable reproduction and high tear resistance are guaranteed.

Triad Custom Tray

For simple fabrication of stable and firm custom impression trays. No waiting required.

Isolant

Dental-stone model and flask isolation liquid for denture fabrication and repairs.

Pinnacle Modelling Wax

Modelling wax perfectly suited for partial and complete denture modelling and tooth set-up and for augmenting bite rims.