



Treatment Centres

Treatment Centres	283
Stools	288
Care and Cleaning Products	289



TREATMENT CENTRES



Treatment Centres



Intego Pro Treatment Centre

Standard Features

For either TS - Traditional Style:
CS - Continental Delivery Systems

- Dentist Element with
 - 1 x 4 Way Sprayvit Syringe - Fibre Optic
 - 1 x Fibre Optic Line
 - 1 x Electric Micro Motor
 - 1 x Sirosonic Piezo Scaler
- Seamless Upholstery
- Articulated Headrest
- Wired Foot Control
- Fully Programmable for 4 users
- Networked & USB Port
- Hand Piece Package
- Dental Assistant Element with
 - 1 x 4 Way Sprayvit Syringe - Fibre Optic
 - 1 x HVE Suction Hose
 - 1 x LVE Suction Hose
- Hand Piece Line Purge and Auto Purge Function
- Manufactured for Wet or Dry Suction Motors and Amalgam Retention
- LED View Plus Operating Light

Optional Extras

- Lounge Upholstery
- Valve Kit - for the connection of External Devices
- Swivelling Tray
- Armrests
- Sivation Package with Camera and Monitor
- Endodontic Function with Reciprocating Endo
- Additional Hand Piece Lines
- 2nd HVE Suction
- Built in Curing Light
- Built in Disinfection and Sanitation for Hand Piece Lines
- Built in Chemical Cleaning for the Suction System
- Heliodent Plus Intra-Oral X-ray Machine
- Ceiling Mounted LED View Plus Operating Light
- Ceiling Mounted LED View Plus Operating Light / Heliodent Plus Combo
- Child Booster Cushion

TENEO

SINIUS

INTEGO PRO

INTEGO



Intego Pro Ambidextrous Treatment Centre

Standard Features

For either TS - Traditional Style : CS - Continental Delivery Systems

- Dentist Element with
 - 1 x 4 Way Sprayvit Syringe - Fibre Optic
 - 1 x Fibre Optic Line
 - 1 x Electric Micro Motor
 - 1 x Sirosonic Piezo Scaler
- Seamless Upholstery
- Articulated Headrest
- Wired Foot Control
- Fully Programmable for 4 users
- Networked & USB Port
- Hand Piece Package
- Dental Assistant Element with
 - 1 x 4 Way Sprayvit Syringe - Fibre Optic
 - 1 x HVE Suction Hose
 - 1 x LVE Suction Hose
- Hand Piece Line Purge and Auto Purge Function
- Manufactured for Wet or Dry Suction Motors and Amalgam Retention
- LED View Plus Operating Light
- Ambidextrous

Optional Extras

- Lounge Upholstery
- Valve Kit - for the connection of External Devices
- Swivelling Tray
- Armrests
- Sivation Package with Camera and Monitor
- Endodontic Function with Reciprocating Endo
- Additional Hand Piece Lines
- 2nd HVE Suction
- Built in Curing Light
- Built in Disinfection and Sanitation for Hand Piece Lines
- Built in Chemical Cleaning for the Suction System
- Heliodent Plus Intra-Oral X-ray Machine
- Ceiling Mounted LED View Plus Operating Light
- Ceiling Mounted LED View Plus Operating Light / Heliodent Plus Combo
- Child Booster Cushion



Intego Treatment Centre

Standard Features

For either TS - Traditional Style : CS - Continental Delivery Systems

- Dentist Element with
 - 1 x 3 Way Sprayvit Syringe or DS Syringe
 - 1 x Fibre Optic Line
 - 1 x Electric Micro Motor
 - 1 x Sirosonic Piezo Scaler
- Seamless Upholstery
- Wired Foot Control
- Flat Headrest
- Fully Programmable for 2 users
- Hand Piece Package
- Dental Assistant Element with
 - 1 x 3 Way Sprayvit Syringe or DS Syringe
 - 1 x HVE Suction Hose
 - 1 x LVE Suction Hose
- Manufactured for Wet or Dry Suction Motors and Amalgam Retention
- LED View Plus Operating Light

Optional Extras

- Lounge Upholstery
- Valve Kit - for the connection of External Devices
- Swivelling Tray
- Armrests
- Articulated Headrest
- Networked & USB Port
- Sivation Package with Camera and Monitor
- Additional Hand Piece Lines
- Suction Hose Cleaning
- Built in Curing Light
- Water Bottle
- Hand Piece Line Purge and Auto Purge Function
- Child Booster Cushion



Sinius Treatment Centre

Standard Features

For either Sliding Track: TS - Traditional Style: or CS - Continental Delivery Systems

- Dentist Element with
 - 1 x 5 Way Sprayvit Syringe - Fibre Optic
 - 1 x Fibre Optic Line
 - 1 x Electric Micro Motor
 - 1 x Sirosonic Piezo Scaler
- Seamless Upholstery
- Articulated Headrest
- Wireless Foot Control
- Fully Programmable for 6 users
- Networked & USB Port
- Hand Piece Package
- Dental Assistant Element with
 - 1 x 5 Way Sprayvit Syringe - Fibre Optic
 - 1 x HVE Suction Hose
 - 1 x LVE Suction Hose
- Built in Disinfection and Sanitation Function for Hand Piece Lines
- Built in Hand Piece Line Purge and Auto Purge Function
- Built in Chemical Cleaning for the Suction System
- Manufactured for Wet or Dry Suction Motors and Amalgam Retention
- LED View Plus Operating Light

Optional Extras

- Lounge Upholstery
- Valve Kit - for the connection of External Devices
- Swivelling Tray
- Armrests
- Motorised Headrest
- Lumbar Support
- Accent Colours
- Sivation Package with Camera and Monitor
- Implant Function
- Endodontic Function with Reciprocating Endo
- Additional Hand Piece Lines
- Surgical Suction
- Built in Curing Light
- Heliodent Plus Intra-Oral X-ray Machine
- Ceiling Mounted LED View Plus Operating Light
- Ceiling Mounted LED View Plus Operating Light / Heliodent Plus Combo
- Child Booster Cushion



Intego Ambidextrous Treatment Centre

Standard Features

For either TS - Traditional Style : CS - Continental Delivery Systems

- Dentist Element with
 - 1 x 3 Way Sprayvit Syringe or DS Syringe
 - 1 x Fibre Optic Line
 - 1 x Electric Micro Motor
 - 1 x Sirosonic Piezo Scaler or Cavitron Scaler
- Seamless Upholstery
- Wired Foot Control
- Flat Headrest
- Fully Programmable for 2 users
- Hand Piece Package
- Dental Assistant Element with
 - 1 x 3 Way Sprayvit Syringe or DS Syringe
 - 1 x HVE Suction Hose
 - 1 x LVE Suction Hose
- Manufactured for Wet or Dry Suction Motors and Amalgam Retention
- LED View Plus Operating Light
- Ambidextrous

Optional Extras

- Lounge Upholstery
- Valve Kit - for the connection of External Devices
- Swivelling Tray
- Armrests
- Articulated Headrest
- Networked & USB Port
- Sivation Package with Camera and Monitor
- Additional Hand Piece Lines
- Suction Hose Cleaning
- Built in Curing Light
- Water Bottle
- Hand Piece Line Purge and Auto Purge Function
- Child Booster Cushion



Hugo Operating Stools

Hugo Range

- 5 Friction Grip Wheels
- 2 x Heights
- Large Range of Upholstery Colours



Hugo Manual
Manual Height Adjustment



Hugo Manual Plus
Manual Height Adjustment with Foot Ring



Hugo Freehand
Height Adjustment by Foot Pressure

Carl Operating Stools

Carl Range

- 5 Friction Grip Wheels
- 2 x Heights
- Large Range of Upholstery Colours



Carl Manual
Manual Height Adjustment
with Belly Bar



Carl Manual Plus
Manual Height Adjustment
with Belly Bar and Foot Ring

Theo Operating Stool

Theo Range

- 5 Friction Grip Wheels
- 2 x Heights
- Large Range of Upholstery Colours



Theo
Manual Height Adjustment with Ring



Dentosept S



Dentosept S (1 L x6) **6610310**
Disinfection Solution of Water Lines for internal use in treatment centres.

Camera Sheaths

UAF Plus Sheath - (box of 300)	6679802
SiroCam C Camera Sheaths (box of 500)	5914705
SiroCam 3, F, AF & AF+ Camera Sheaths (box of 500)	5938654

Quality you can rely on

We don't want to just satisfy our customers, we want to impress them. This is why quality, performance, and service are key for us. Dentsply Sirona treatment units combine high quality with real innovations that perceptibly enhance treatment comfort and workflow.

Reliability and durability



Precise, elaborate quality tests and robust materials make Intego Pro a reliable treatment center with an above-average service life.

High quality



Intego Pro immediately shows our attention to detail. The materials used – robust, shapeable plastics, solid steel and aluminum components on the inside, high-quality glass and upholstery – meet the high quality demands of a premium product.

High standards



Dentsply Sirona is a technology leader with more than 130 years of innovation history. Our engineers and scientists in Research & Development not only pay close attention to innovation, but also quality.



Future-proof design

At the Dentsply Sirona innovation center, future trends are anticipated and products are prepared to be equipped for the technology demands of the future.

